



Round Handle Cane

The **ProBasics Round Handle Cane** provides excellent support and is designed for patients who need additional assistance with balance. The round handle allows users to hook the cane over their arm while opening doors or performing other tasks. Supports patient weights up to 250 pounds.



Easy one-button height adjustment and anti-rattle ring



Round vinyl handle

FEATURES

- Easy one-button height adjustment with a locking anti-rattle ring allowing for quick, safe adjustment
- Adjusts in 1-inch increments
- Anodized aluminum finish with vinyl grip
- Curved handle allows cane to be retained on arm
- Easy-to-clean aluminum frame
- HCPCS* Code: E0100

SPECIFICATIONS

Handle Height Adjustment Range:	31" – 39"
Cane Diameter:	7/8"
Cane Tip:	3/4"
Weight Capacity:	250 lb
Product Weight:	13.44 oz
Carton Dimensions:	32" x 6.7" x 8.0"
Shipping Weight:	9.24 lb
Warranty:	Limited Lifetime on frame 90 days on all parts

Item #	Description	Color	UPC	UOM
CNRSL	ProBasics Round Handle Cane	Silver	815067071265	10/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.