



## QUICK TIPS

- This device has been thoroughly tested prior to delivery including heating performance verification and leak tests.
- Distilled water helps to maintain device cleanliness. The user guide recommends 12 gallons of distilled water and 1 cup of tap water.
- Discoloration of the temperature sensor is expected and results from the normal interaction between the electroless-plated nickel surface and elevated water temperatures.
- For routine cleaning and maintenance, use a soft lint free cloth with isopropyl alcohol. Do not use abrasive pads or harsh chemicals. Refer to the User Guide for additional guidance on cleaning routines.