



Carex OTC Finger Pulse Oximeter

The **Carex OTC Finger Pulse Oximeter** allows you to monitor and maintain pulse and oxygen saturation levels for optimal health. Accurate, fast oxygen level and heart rate measuring technology is ideal for sports enthusiasts. Recommended for spot-check monitoring only. No prescription required.



FEATURES

- Approved for over-the-counter retail sale
- Large easy-to-read display
- Suitable for pediatric and adult use
- Soft silicone padding for a comfortable fit
- One-button operation
- Up to 30 hours of use with batteries
- Auto-off function
- Low battery indicator
- Manual, lanyard and two AAA alkaline batteries included
- Attractive retail packaging
- Two-year Limited Warranty

SPECIFICATIONS

Product Dimensions:	2.25" L x 1.25" W x 1.5" H
Product Weight: <i>(including two AAA alkaline batteries)</i>	1.8 oz
Display:	LED
Display Range:	0-99% (SpO ₂) / 0~254 BPM
Measurement Range:	70-99% (SpO ₂) / 30-235 BPM
Accuracy:	70-99%: ±3% (SpO ₂) / 30~99 BPM: ±2BPM; 100~235 BPM: ±2%
Resolution:	±1% for SpO ₂ / ±1 BPM for Pulse Rate
Operation Temperature:	41~105°F (5~40°C)
Storage Temperature:	-4~131°F (-20~55°C)
Ambient Humidity:	15%~80% (operation); ≤93% (storage)

Item #	Description	UPC	UOM
FGR41800 0000	Carex OTC Finger Pulse Oximeter	023601064189	each