



## Standard Rollator with 6-inch Wheels

The **ProBasics Aluminum Rollator** is designed to fit the majority of your patients' fundamental needs with a focus on style and affordability. Constructed of lightweight aluminum, it offers great maneuverability and stability. Equipped with easy-to-operate ergonomic hand brakes, a removable padded backrest and a water-resistant nylon pouch. Easily folds for storage and transport with built-in fold strap. Supports patient weights up to 300 pounds.



Water-Resistant  
Nylon Pouch



Easy-to-Operate  
Ergonomic Hand Brakes



6-inch Wheels

### FEATURES

- Lightweight aluminum frame in burgundy, blue, and black
- Two-button removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 6" black casters with soft grip tires, ideal for both indoor and outdoor use
- HCPCS\* Code: E0143  
E0156 (seat)

### COLORS

Burgundy Finish  Blue Finish  Black Finish 

### SPECIFICATIONS

Overall Width:	24"
Width b/w Hand Grips:	19"
Seat Dimensions:	14" (W) x 12" (D)
Seat-to-Floor Height:	21"
Weight Capacity:	300 lb
Adjustable Handle Height:	30" – 36"
Wheel Size:	6"
Product Weight:	15.4 lb
Carton Dimensions:	25" L x 8" W x 34" H
Shipping Weight:	16 lb
Warranty:	Limited lifetime on frame

Item #	Description	Color	UPC	UOM
RLA6BG	ProBasics Standard Rollator with 6-inch Wheels	Burgundy	815067071630	1/cs
RLA6BL	ProBasics Standard Rollator with 6-inch Wheels	Blue	815067071647	1/cs
RLA6BK	ProBasics Standard Rollator with 6-inch Wheels	Black	815067071654	1/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.