



# Aspen Dual Height Adjustable Rollator

Experience comfort, stability, and convenience with the Aspen Dual Height Adjustable Rollator. Crafted from lightweight aluminum, oval-shaped tubing and 8" wheels for both indoor and outdoor use, the Aspen is both attractive and built for years of worry-free use. This dual adjustable rollator features four seat height and five handle height adjustments, making it easy to achieve the perfect ergonomic fit for a wide variety of users. Comfort is emphasized with the Aspen's large, comfortable padded seat and extra-padded strap-style backrest that is designed to support the user with superior comfort through rest breaks. With first-time tool-free assembly, setup is quick and hassle-free. See why the Aspen is the best choice for both the dealer and the user for Rollator design, comfort, durability, and value.

## SPECIFICATIONS

<b>Weight Capacity:</b>	<b>300 lb</b>
Overall Width:	23"
Seat Dimensions:	14.5" (W) x 14.5" (D)
Patient Height Range:	5'2" - 6'3"
Handle Height Range:	
<b>Seat Height Setting</b>	<b>Handle Height Range</b>
Seat-to-Floor Height: 20"	34" - 38"
Seat-to-Floor Height: 21"	35" - 39"
Seat-to-Floor Height: 22"	36" - 40"
Seat-to-Floor Height: 23"	37" - 41"
Wheel Size:	8"
Product Weight:	16.4 lbs
Carton Dimensions:	20" x 8" x 22"
Shipping Weight:	17 lbs
Limited Lifetime Warranty:	5 Year Warranty on Frame 1 Year on Parts

## FEATURES

- Features an attractive, lightweight aluminum oval-tube frame that folds easily for storage and on the go travel – a perfect design for years of worry-free use
- 8" all terrain wheels are perfect for both indoor and outdoor use
- (4) 1" Seat height and (5) 1" handle height adjustments to custom fit a variety of users
- Extra padded flip-up seat offers excellent comfort and easy access to an under-seat storage bag
- New flip-up backrest that is flexible, wider, more comfortable, and height adjustable
- New easy-adjustable rear braking system with covered calipers
- Tool Free first-time assembly is fast and simple
- Weight capacity: 300 lbs
- Warranty: Limited lifetime on frame
- HCPCS\* Code: E0143, E0156 (seat)

## COLORS



WARNING: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer or birth defects or other reproductive harm. For more information, visit [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Item #	Description	UPC
RLAH08CB	Aspen Rollator, Cobalt Blue	815067012695
RLAH08BW	Aspen Rollator, Brandywine	815067012701
RLAH08BK	Aspen Rollator, Black	815067012725

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.