



Ascent Euro-Style Rollator

The Ascent is not only our newest model, it's the Flagship of our robust line of rollators. This lightweight, euro-style, 4-wheel rollator has both an adjustable height padded seat and backrest. The all new Ascent boasts a number of high-end features to ensure safety and longevity including dual sided, vented aluminum front forks, stylish 10" front and 8" rear casters with black non-marring tires and a unique top-to-bottom protected braking system. The Ascent is featured in three brilliant colors: Cobalt Blue, Brandywine, and Charcoal Black finish.



- Easy push button for handle height adjustments in 1" increments
- Aluminum brake cable conduit protects at the most frequent breakage point



- Innovative design hides brake lines for a clean look while protecting from wear and breakage
- Shielded adjustable braking mechanism offers excellent durability and performance




- Dual Sided Aluminum Front Forks with reinforced 10" Front & 8" rear casters for longevity and durability

FEATURES

- Lightweight Aluminum Frame for ease of use and longevity.
- Ergonomic handles with hidden brake cables; adjustable height in one-inch increments from 37.5" – 42.5" with an audible click as well as visual push button height settings.
- New and improved Dual Sided Aluminum Front Forks for increased stability.
- 18" x 10" thickly padded seat with a handle to quickly and easily fold and unfold side-to-side (rather than front-to back), allowing the unit to maneuver through narrow doorways.
- Dual Sided, vented aluminum front forks with 10-inch front and 8-inch rear heavy-duty wheels with black non-marring, soft-grip tires to traverse almost any terrain.
- Top to Bottom hidden braking system is fully adjustable, durable and stylish.

COLORS



SPECIFICATIONS

Weight Capacity:	300 lb
Overall Width:	23.5" (W)
Seat Dimensions:	18" (W) x 10" (D)
Seat Height:	Handle Height Range:
21"	37.5" - 40.5"
22"	38.5" - 41.5"
23"	39.5" - 42.5"
Patient Height Range:	5'6" - 6'5"
Adjustable Handle Height:	37.5" – 42.5"
Wheel Size:	10" (front) 8" (rear)
Product Weight:	16.4 lb
Carton Dimensions:	20.7" x 8.85" x 26"
Shipping Weight:	24.5 lb
Limited Lifetime Warranty:	5 Year Warranty on Frame 1 Year on Wear Parts

* User body types vary, heights are approximate based on average person.

WARNING: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer or birth defects or other reproductive harm. For more information, visit www.P65Warnings.ca.gov

Item #	Description
RLAS10CB	Ascent Rollator, Cobalt Blue
RLAS10BW	Ascent Rollator, Brandywine
RLAS10BK	Ascent Rollator, Textured Black

* Map Pricing: \$259 Cannot be sold on Amazon

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.