




# Voyager™ Euro-Style Rollator

The Voyager is a unique US/European hybrid design with separate seat and handle height ranges that accommodate heights ranging from 4'11" – 6'4" and as the seat height is increased the footprint grows effectively addressing the need for increased stability as the center of gravity of the unit moves upward. This lightweight rollator features easy-to-operate ergonomic hand brakes with hidden cables. It has 10-inch front and 8-inch rear wheels for great maneuverability and stability with an attractive, shielded braking mechanism which offers excellent grip on the rear wheels.



- Easy push button for handle height adjustments in 1" increments
- Aluminum brake cable conduit protects at the most frequent breakage point



- Innovative design hides brake lines for a clean look while protecting from wear and breakage
- Shielded adjustable braking mechanism offers excellent durability and performance



- 10-inch Front Caster Reduces Resistance for Easy Navigation

## FEATURES

- Ergonomic handles with hidden brake cables; adjustable height in one-inch increments from 35" – 42.5" with an audible click as well as visual height settings.
- Height adjustments on leg frame allows seat height to be adjusted from 21" to 24".
- Comes with an 18" x 10" thickly padded seat with a handle to quickly and easily fold and unfold side-to-side (rather than front-to-back), allowing the unit to maneuver through narrow doorways.
- Attractive mesh tote attaches to the front frame and is easily removed to double as a shoulder bag for carrying personal items.
- 10-inch front and 8-inch rear non-marring black wheel with soft grip tires glide on almost any surface.
- Lightweight aluminum frame is durable for years of worry-free use.
- Folds up for quick and easy transport in any vehicle and for storage.

## SPECIFICATIONS

Setting	Low	Mid	High
Seat to Floor Height	21"	22"	24"
Overall Height Range	35"– 39.5"	36"– 40.5"	38"– 42.5"
<i>Brake Handles adjust between 14"– 18.5" from the seat, in 1" increments</i>			
Overall Width:	23.5" (W)		
Seat Dimensions:	18" (W) x 10" (D)		
Adjustable Seat Height:	21"– 24"		
Patient Height Range:	4'11" - 6'4"		
Weight Capacity:	300 lb		
Adjustable Handle Height:	35"– 42.5"		
Wheel Size:	10" (front) 8" (rear)		
Product Weight:	16.4 lb		
Carton Dimensions:	26" x 8.85" x 20.7"		
Shipping Weight:	21 lb		
Warranty:	Limited lifetime on frame		

## COLORS

- Cobalt Blue  Rose Gold  Ice Palace 

Item #	Description	Color	UPC	UOM
RLEU10BL	Voyager Euro-Style Rollator	Cobalt Blue	815067078066	1/cs
RLEU10PK	Voyager Euro-Style Rollator	Rose Gold	815067078073	1/cs
RLEU10WT	Voyager Euro-Style Rollator	Ice Palace	815067078080	1/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.