



Knocked Down Economy Steel Rollator with 6-inch Wheels

The **ProBasics Steel Rollator** is the perfect blend of style, comfort and affordability. Constructed of durable powder-coated steel, it offers great maneuverability and stability. It features easy-to-operate ergonomic hand brakes, a removable padded backrest and a water-resistant nylon pouch. Equipped with an easy folding mechanism for storage and transport. Supports patient weights up to 300 pounds. The new knock-down version features a quick tool-free first time assembly.



Water-Resistant Nylon Pouch



Easy-to-Operate Ergonomic Hand Brakes



6-inch Wheels

FEATURES

- Durable powder-coated steel frame
- Removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 6" black casters with soft grip tires, ideal for both indoor and outdoor use
- Easy Tool-Free setup - Installation guide included.
- New & Improved Box Size for Increased Product Availability reduced shipping costs.
- HCPCS* Code: E0143

COLORS

Burgundy
Finish



Blue
Finish



SPECIFICATIONS

Overall Width:	23"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13" (W) x 9.5" (D)
Seat-to-Floor Height:	20.5"
Weight Capacity:	300 lb
Adjustable Handle Height:	30.5" - 36.5"
Wheel Size:	6"
Product Weight:	16.5 lb
Carton Dimensions:	23.0" x 4.75" x 22"
Shipping Weight:	19 lbs
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLS6BG-KD	ProBasics Economy Steel Rollator with 6-inch Wheels	Burgundy	815067079926	1/cs
RLS6BL-KD	ProBasics Economy Steel Rollator with 6-inch Wheels	Blue	815067079919	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.