

Voyager XR[™] **Euro-Style Rollator**

The Voyager XR is a unique extra-rugged rollator that accommodates a wide range of users. One of the differentiating features of the rollator is that has a larger weight compacity of 450 lbs. It also features a rear shock absorbing system, which greatly improves our user's stability and comfort as they operate the rollator as well as both sit in and stand up from sitting in the seat. One of the best features of the Voyager XR is the top-to bottom protected braking system. With a shielded braking mechanism, ergonomic hand grips, and hidden cable management, the Voyager XR truly has it all!



handle height adjustments

conduit protects at the most

frequent breakage point

in 1" increments

• Aluminum brake cable



 Innovative design hides brake lines for a clean look while protecting from wear and breakage

 Shielded adjustable braking mechanism offers excellent durability and performance



- 10-inch Front Caste Reduces Resistance for Easy Navigation

SPECIFICATIONS

Overall Width:	29" (W)		
Seat Dimensions:	24" (W) x 13" (D)		
Seat to Floor Height:	24″		
Weight Capacity:	450 lb		
Adjustable Handle Height:	36"-41"		
Wheel Size:	10" (front) 8" (rear)		
Product Weight:	23 lb		
Carton Dimensions:	29.7" x 11.2" x 37.6"		
Shipping Weight:	30.4 lb		
Warranty:	Limited Lifetime on Frame Wear Components: 1 year		

FEATURES

- Lightweight, aluminum frame with integrated suspension system is durable for years of worry-free use.
- 10-inch front and 8-inch rear non-marring, heavy duty black wheels with soft grip tires traverse almost any terrain.
- Ergonomic handles with top-down hidden braking system; Adjustable height in one-inch increments from 36"- 41" with audible "click" setting.
- Comes with a 24" x 13" thickly padded seat with a handle to fold and unfold side-to-side quickly and easily. This allows for the rollator to maneuver through narrow hallways and doors.
- Attractive front mesh tote attaches to the front of the frame for carrying personal belongings. It is easily attached and detached and can double as a shoulder bag when necessary.
- Folds up quickly and easily for transport or storage. Can easily fit in vehicles, storage closet, or public transport.

COLORS



ltem #	Description	Color	UPC	UOM
RLVXR10CB	Voyager XR Euro-Style Rollator	Cobalt Blue	815067078646	1/cs
RLVXR10RG	Voyager XR Euro-Style Rollator	Rose Gold	815067078639	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.