



Rite-Neb 5

NEBULIZER COMPRESSOR SYSTEM

The Roscoe RITE-NEB5 Nebulizer System offers reliability, high performance and convenient size with a convenient carry handle. The simple one-button operation and lightweight design makes it extremely portable and patient friendly. Weighing under 3 lbs, the Rite-Neb 4 includes the compressor, nebulizer, tee, mouthpiece, corrugated reservoir tube, tubing, five extra filters, instruction manual and a five year warranty. The Roscoe RITE-NEB5 is also available with a reusable neb kit option.

PRODUCT FEATURES:

- One-button operation
- Delivers better output for faster treatments and better patient compliance
- Compact design makes it portable
- Quiet at less than 65 dBA
- Compact and lightweight at 2.7 lbs
- Metal piston for durability
- Attractive retail packaging
- Comes complete with compressor, nebulizer, tee, mouthpiece, reservoir tube, tubing and five extra filters
- TruNeb Reusable Neb Kits available – purchase with the compressor or alone
- HCPCS Code: E0570*
- HCPCS Code: A7003* for nebulizer kits
- HCPCS Code: A7005* for reusable nebulizer kits
- HCPCS Code: A7013* for felt filters

PRODUCT SPECIFICATIONS:

Dimensions:	6.1" L x 5.5" W x 3.5" H
Weight:	2.7 lbs
Method of Nebulization:	Piston compressor
Maximum Pressure:	30 to 50 psi
Operating Pressure Range:	9.5 PSI - 16 PSI
Maximum Flow:	6 - 10 LPM
Sound Level:	<65 dBA
Medication Capacity:	5 mL
Particle Size:	< 0.5 µm
Electrical Requirements:	120 VAC/60Hz
Operation Cycle:	20 min. on, 30 min. off
Warranty:	Limited 5 year

OPTIONAL ACCESSORIES:



Item #	Description	Qty
RITE-NEB5	RITE-NEB5 Nebulizer Compressor System with Disposable Neb Kit	each
RITE-NEB5-REUS	RITE-NEB5 Nebulizer Compressor System with Reusable Neb Kit	each
RITE-NEB5-REUSPM-PB	RITE-NEB5 Nebulizer Compressor System with one disposable neb kit, one reusable neb kit, one pediatric mask with piggyback serial labels	each
Optional Accessories:		
CIF-1128PK6	Foam Filter, 5/8" diameter x 1/2" thick	6/pk
50016	TruNeb Reusable Nebulizer Kit	10/cs
BAG-NEBR	Nylon Carry Bag	each
BAG-DRW	Nylon Drawstring Carry Bag	each

* This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional. Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.