

## ASCENT XR ADJUSTABLE HEIGHT EURO-STYLE ROLLATOR

RLASXRCB/RG/BW



**COLORS** Cobalt Blue Brandywine Finish

- Lightweight, aluminum frame with integrated suspension system is durable for years of worry-free use.
- Ergonomic handles – with top-down hidden braking system; Adjustable height in one-inch increments from 36" - 41" with audible "click" setting.
- Includes a densely padded adjustable height backrest for maximum comfort and support.

### SPECIFICATIONS

Overall Width:	29" (W)
Seat Dimensions:	24" (W) x 13" (D)
Seat to Floor Height:	24"
Weight Capacity:	450 lb
Adjustable Handle Height:	36" – 41"
Wheel Size:	10" (front) 8" (rear)
Product Weight:	23 lb
Carton Dimensions:	29.7" x 11.2" x 37.6"
Shipping Weight:	30.4 lb

Warranty: Limited lifetime on frame. One year on all parts.

## ASCENT ADJUSTABLE HEIGHT EURO-STYLE ROLLATOR

RLAS10CB/BK/BW



**COLORS** Cobalt Blue Black Finish Brandywine Finish

- Includes 10-inch front and 8-inch rear wheels.
- New and improved Dual Sided Aluminum Front Forks for increased stability.
- Includes a densely padded adjustable height backrest for maximum comfort and support.
- Top to Bottom hidden braking system is fully adjustable, durable and stylish.

### SPECIFICATIONS

Overall Width:	23.5" (W)
Seat Dimensions:	18" (W) x 10" (D)
Adjustable Seat Height:	21" – 24"
Weight Capacity:	300 lb
Adjustable Handle Height:	37.5" – 42.5"
Wheel Size:	10" (front) 8" (rear)
Product Weight:	16.4 lb
Carton Dimensions:	20.7" x 8.85" x 26"
Shipping Weight:	17.5 lb

Warranty: Limited lifetime on frame. One year on all parts.

## ASPEN ADJUSTABLE HEIGHT ALUMINUM ROLLATOR

RLAH08BL/BK/BW/RG



**COLORS\*** Blue Finish Black Finish Brandywine Finish Rose Gold

- First-time tool-free assembly required.
- Extra-padded strap style back for maximum comfort.
- 5 - 1" Height adjustments to fit a wide range of users.
- 4 - 1" Seat height adjustments for improved seat-to-floor height.
- Featured in Blue, Black, Brandywine, and Rose Gold\*.

### SPECIFICATIONS

Overall Width:	23" (W)
Seat Dimensions:	14" (W) x 12" (D)
Adjustable Seat Height:	18" – 22"
Weight Capacity:	300 lb
Adjustable Handle Height:	30" – 38"
Wheel Size:	8"
Product Weight:	14 lb
Carton Dimensions:	31.5" x 7.5" x 24.6"
Shipping Weight:	19 lb

\*Colors subject to change

## KNOCKED-DOWN ECONOMY STEEL ROLLATOR WITH 6-INCH WHEELS

RLS6BG-KD, RLS6BL-KD



**COLORS** Burgundy Finish Blue Finish

- Requires a 5-minute, tool-free first-time assembly for increased cost savings.
- Includes non-marring 6" black caster wheels with soft grip tires.
- Includes a 13" x 9.5" padded seat and thick padded backrest for increased comfort.
- Weight capacity of up to 300 lbs.
- Featured in Burgundy and Blue.

### SPECIFICATIONS

Overall Width:	23"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13" (W) x 9.5" (D)
Seat-to-Floor Height:	20.5"
Weight Capacity:	300 lb
Adjustable Handle Height:	30.5" – 36.5"
Wheel Size:	6"
Product Weight:	16.5 lb
Carton Dimensions:	23" x 4.75" x 22"
Shipping Weight:	19 lb

Warranty: Limited lifetime on frame. One year on all parts.

## PROBASICS ALUMINUM ROLLATOR WITH 6-INCH WHEELS

RLA6BG/BL/BK



**COLORS** Burgundy Finish Blue Finish Black Finish

- A lightweight frame that includes non-marring 6" black casters with soft grip tires, ideal for both indoor and outdoor use.
- Weight capacity of up to 300 pounds.
- Featured in Blue, Black, and Burgundy.

### SPECIFICATIONS

Overall Width:	23"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13" (W) x 9.5" (D)
Seat-to-Floor Height:	20.5"
Weight Capacity:	300 lb
Adjustable Handle Height:	30.5" – 36.5"
Wheel Size:	6"
Product Weight:	16.5 lb
Carton Dimensions:	24.6" x 7.5" x 31.5"
Shipping Weight:	19 lb
Warranty: Limited lifetime on frame. One year on all parts.	

## PROBASICS DELUXE ALUMINUM ROLLATOR WITH 8-INCH WHEELS

RLA8BL/BG/BK



**COLORS** Burgundy Finish Blue Finish Black Finish

- A lightweight frame that is equipped with large 8-inch wheels for improved maneuverability on outdoor terrain.
- Weight capacity of up to 300 pounds.
- Featured in Blue, Black, and Burgundy.

### SPECIFICATIONS

Overall Width:	22.75"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13.75" (W) x 11.5" (D)
Seat-to-Floor Height:	21"
Weight Capacity:	300 lb
Adjustable Handle Height:	31" – 36"
Wheel Size:	8"
Product Weight:	15.4 lb
Carton Dimensions:	24.8" x 9.1" x 35"
Shipping Weight:	18.3 lb
Warranty: Limited lifetime on frame. One year on all parts.	

## BARIATRIC ROLLATOR WITH 8-INCH WHEELS

RLAB8BG/BL



**COLORS** Burgundy Finish Blue Finish

- Made with an extra-strength, reinforced aluminum frame.
- Includes non-marring 8" black casters with reinforced forks and soft grip tires, ideal for both indoor and outdoor use.
- Weight capacity of up to 400 pounds.
- Featured in Burgundy, and Blue.

### SPECIFICATIONS

Overall Width:	29"
Width b/w Hand Grips:	23.5"
Seat Dimensions:	19" (W) x 13" (D)
Seat-to-Floor Height:	22"
Weight Capacity:	400 lb
Adjustable Handle Height:	33.5" – 38.5"
Wheel Size:	8"
Product Weight:	18.7 lb
Carton Dimensions:	29.1" x 8.6" x 36.2"
Shipping Weight:	22.7 lb
Warranty: Limited lifetime on frame. One year on all parts.	

## PROBASICS JUNIOR ROLLATOR WITH 6-INCH WHEELS

RLAJ6BG/BL



**COLORS** Burgundy Finish Blue Finish

- Designed to fit individuals of a shorter stature.
- Two-button removable curved padded backrest.
- Comfortable, flip-up padded seat.
- Easy-to-operate ergonomic hand brakes.
- Weight capacity of up to 250 pounds.
- Featured in Burgundy, and Blue.

### SPECIFICATIONS

Overall Width:	22"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13" (W) x 13.5" (D)
Seat-to-Floor Height:	19"
Weight Capacity:	250 lb
Adjustable Handle Height:	29" – 34"
Wheel Size:	6"
Product Weight:	12 lb
Carton Dimensions:	22.6" x 7.3" x 31.3"
Shipping Weight:	14.3 lb
Warranty: Limited lifetime on frame. One year on all parts.	