

Product Category Description: Oxygen Concentrator	Item #
truAIRE-5 Stationary 5-Liter Oxygen Concentrator	O2C5L

	Replacement Parts for: <b>O2C5L</b>	CHB Item#
1	OXYGEN OUTLET FITTING	X1JAO2C01
2	FLOW METER WITH KNOB	X1JAO2C02
3	HOUR METER	X1JAO2C04
4	POWER SWITCH	X1JAO2C05
5	CIRCUIT BREAKER	X1JAO2C06
6	CASTER, 4/PK	X1JAO2C07
7	POWER CORD	X1JAO2C08
8	CABINET FILTER PLATE <i>(not pictured)</i>	X1JAO2C09
9	SIEVE BEDS, PAIR	X1JAO2C10
10	CONTROL VALVE	X1JAO2C11
11	PRESSURE REGULATOR	X1JAO2C12
12	PE VALVE	X1JAO2C13
13	PRODUCT TANK	X1JAO2C14
14	UNIT BASE	X1JAO2C15
15	OXYGEN SENSOR	X1JAO2C16
16	CAPACITOR	X1JAO2C17
17	PC BOARD	X1JAO2C18
18	TRANSFORMER	X1JAO2C19
19	EXHAUST FAN	X1JAO2C20
20	DOUBLE WOBBLE COMPRESSOR	X1JAO2C21
21	COMPRESSOR BRACKET	X1JAO2C22
22	CONTROL PANEL	X1JAO2C23
23	COMPRESSOR CUP SEAL	X1JAO2C24
24	COMPRESSOR REED VALVE	X1JAO2C25
25	COMPRESSOR GASKET	X1JAO2C26
26	PRODUCT TANK CHECK VALVE	X1JAO2C27
27	EXHAUST MUFFLER	X1JAO2C28



	Replacement Parts for: <b>O2C5L</b>	CHB Item#
28	INTAKE RESONATOR	X1JAO2C29
29	WIRE HARNESS (O2 SENSOR)	X1JAO2C30
30	WIRE HARNESS (MAIN POWER)	X1JAO2C31
31/32	REPLACEMENT EXTERIOR CABINET (31-FRONT, 32-BACK)	X1JAO2C33
33	SHIPPING CARTON ( <i>not pictured</i> )	X1JAO2C34
34	BASE FOAM INSERT ( <i>not pictured</i> )	X1JAO2C35
35	TOP FOAM INSERT ( <i>not pictured</i> )	X1JAO2C36
36	CONCENTRATOR REPAIR LABOR ( <i>not pictured</i> )	X1JAO2C37
37	EXHAUST AIR GUIDE COVER	X1JAO2C38
38	MISC CONC PART ( <i>not pictured</i> )	X1JAO2C39
39	CABINET FILTER, 20/PK	CF-1000
40	INTAKE FILTER	CIF-O2C5L
41	FINAL BACTERIA FILTER	ILB-316



