



Meridian Medical Air-Foam APM System



The **Meridian SatinAir Foam (5" air/3" foam) APM system** is a complete therapy mattress system. Designed to provide alternating pressure and low air loss, the SatinAir Foam optimizes pressure redistribution and helps to manage skin maceration. This system is the perfect option where a base foam layer is critical to prevent patients bottoming out during power outages.

The powerful 9 liter-per-minute pump offers static inflation for transfers or patient care and convenient settings by patient weight. The nylon zipper cover is low shear, fluid-resistant and vapor permeable. The emergency CPR valve allows for quick deflation if necessary.

The **Meridian SatinAir Foam** offers innovative design for use in nursing homes, hospitals and the home care environment.

FEATURES

- Includes 5" air cell/3" foam mattress, cover and pump
- Pillow function maintains the air in cells at the head for patient comfort
- Air cells can be easily removed and replaced
- CPR quick deflate at head of mattress
- Quilted nylon zipper cover is low shear, fluid-resistant, and vapor permeable
- 9 liter-per-minute pump in an easy-to-use design
- Two modes - static or alternating
- Variable pressure setting for maximum comfort

MATTRESS CONSTRUCTION

- Seventeen rows of 5" alternating pressure air cells, 9 with laser holes
- 3" foam base to prevent bottoming out during power outages
- Air cells of durable Nylon/PVC for easy cleaning or removal and replacement

COVER

- Zippered for convenience
- Quilted, vapor-permeable nylon for patient comfort
- Fluid-resistant
- Easy cleaning
- Low shear/low friction

PRODUCT SPECIFICATIONS

Model:	SA-SAF350A
Mattress:	36" wide x 8" high x 80" long
Pump:	9 LPM analog, 2 modes
Cells:	Seventeen 5" cells, nylon/PVC with 3" foam base
Patient Weight:	350 lbs. maximum
Pressure Range:	30 mmHg – 60 mmHg (+/- 10 mmHg)
Product Weight:	30.9 lbs.
Shipping Dimensions:	18" x 38.5" x 20", 33 lbs.
Warranty:	Limited, 1 year on pump and mattress
HCPS Code:	E0277*

PUMP

- 9 liter-per-minute analog pump
- Front panel switch provides two modes - static or alternating
- Patient weight setting on the dial for quick set up
- Audible low pressure alarm

Item #	Description	UOM
SA-SAF350A	Meridian Medical Air-Foam APM System	each

* - Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is only included for your convenience and is not intended, nor should be considered as Medicare billing advice.