

# Sliding Transfer Bench with Swivel Seat



The **Sliding Transfer Bench with Swivel Seat** provides an easy way to transfer in and out of the tub. The swivel seat turns 360° and locks every 90° allowing for easy entry and exit from the bench. The seat slides smoothly left and right along the frame, and comes equipped with a seat belt for added safety.



Seat swivels and slides left to right for easy transfers



Versatile to accommodate any bathroom

## FEATURES

- Swivel seat turns 360° and locks every 90° allowing easy entry to and exit from the bench
- Specially-designed sliding rings lock the seat at each end of the sliding rails
- Seat glides over polished, high-strength aluminum tubes
- Molded plastic seat and back, textured finish
- Rust-proof, lightweight aluminum construction
- Adjustable height to fit user and bathtub clearance
- Comes with safety belt strap

## SPECIFICATIONS

Width (Outside Legs):	39.5"
Back Height (Outside Legs):	13.75"
Seat Width:	17.75"
Seat Depth:	14.75"
Adjustable Seat Height:	18.75" - 23"
Weight Capacity:	400 lb
Warranty:	One year manufacturer warranty

Item #	Description	UPC	UOM
77662	Sliding Transfer Bench with Swivel Seat	604180776627	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.