

Sliding Transfer Bench with Seat and Back

The Sliding Transfer Bench with Seat and Back provides an easy way to transfer in and out of the tub. The bench slides smoothly left and right along the frame, and comes equipped with a seat belt for added safety.





Seat slides left to right for easy transfers

Versatile to accomodate

FEATURES

- Specially-designed sliding rings lock the seat at each end of the sliding rails
- Seat glides smoothly in and out of the bathtub or shower over polished, high-strength aluminum tubes
- Comfortable, molded plastic seat with textured finish
- Rust-proof, lightweight aluminum construction
- Adjustable height to fit user and bathtub clearance
- Comes with safety belt strap
- No tools necessary for assembly

SPECIFICATIONS

Width (Outside Legs):	39.5"	
Back Height (Outside Legs):	13.75"	
Seat Width:	17.75"	
Seat Depth:	14.75"	
Adjustable Seat Height:	18.25" - 23"	
Weight Capacity:	400 lb	
Warranty:	One year manufacturer warranty	

Item #	Description	UPC	UOM
77211	Sliding Transfer Bench with Seat and Back	604180772117	each

^{*} Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.