

TheraSonic™

CONDUCTIVE GEL

WATER-SOLUBLE CONTACT MEDIUM FOR ULTRASOUND TRANSMISSION

- Acoustically correct to accurately transmit the broad range of ultrasound frequencies
- Hypoallergenic, bacteriostatic & non-irritating
- Formaldehyde-free

DIRECTIONS FOR USE

- Apply gel to treatment area
- Conduct treatment
- Remove with tissue or towel
- Clean equipment after each use

FOR EXTERNAL USE ONLY.
NOT FOR USE WITH DEFIBRILLATORS.



Manufactured for:
Compass Health Brands Corp.
4120 South Creek Road
Chattanooga, TN 37406



Richmar®

3120.18.01A ©2018 Compass Health Brands Corp.

Ingredients:

Water, EDTA Disodium,
Propylene Glycol,
Triethanolamine,
Carbomer, FD&C Blue
#1, Diazolidinyl Urea,
Iodopropynyl
Butylcarbamate

Made in Canada

16 FL OZ
(132 USG / 473 mL)

EXP. DATE:
YYYY-MM-DD

LOT#: