



K3 Lite Wheelchair

LIGHTWEIGHT WHEELCHAIR

The **ProBasics K3 Lite Wheelchair** is a standard lightweight, extremely versatile wheelchair. The K3 Lite features a stylish, easy-to-clean nylon upholstery and comes standard with a variety of different features and benefits. The K3 Lite easily adapts to any patient's individual needs. It offers superior value and functionality for those users looking for a high-performance mobility product at an affordable price.



Flip-back armrests are height adjustable



2" upholstery seat extension



Adjustable anti-tippers with push-pin flip-up wheels



24" Rear Wheels

FEATURES

- High-Strength Steel Frame for years of worry free use
- Comes included with a seat rail and upholstery extension that easily adjusts seat depth from 16" to 18" without any additional upholstery
- Flip-back padded arm rests are height adjustable from 10" to 14"
- 3-Position axles allow for removal of rear wheels for fast hemi-height adjustment
- Comfortable, lightweight and stylish easy-to-clean nylon upholstery
- 24" Rear Polyurethane tires and 8" composite casters for superior performance
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Includes tool-free adjustable front riggings
- Easy push-to-lock wheel locks - Brake Extensions Included in Legrest Box.
- HCPCS* Code: K0003
K0195 or E0990 (elevating legrests only)
E0971 (anti-tippers)

SPECIFICATIONS

Seat Width:	16" 18" 20"
Seat Depth:	16" – 18"
Seat-to-Floor Height ¹ :	Adult 19.5" Hemi 17.5"
Back Height:	17"
Back Style:	Fixed
Overall Width w/Arms:	24.5" 26.5" 28.5"
Overall Depth:	32" (without riggings)
Weight Capacity:	300 lb
Product Weight ² :	35 - 38 lbs
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	49 – 55 lbs (Includes Legrests)
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC3L1616DS	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067012763	each
WC3L1816DS	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067012787	each
WC3L2016DS	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067012800	each
WC3L1616DE	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067012770	each
WC3L1816DE	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067012794	each
WC3L2016DE	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067012817	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.