



K4 Deluxe Wheelchair

HIGH STRENGTH K4 WHEELCHAIR

The **ProBasics K4 Deluxe Wheelchair** is a high-strength, feature rich wheelchair. This wheelchair comes standard with a 300 LB weight capacity, One -Button-Convertible Desk-To-Full Length Armrests, Quick Release Axle for Hemi-Height Adjustments. The K4 Deluxe easily adapts to any patient's individual needs. It offers superior value and functionality for those looking for a high-performance mobility product.



2" upholstery seat extension



Adjustable anti-tippers with push-pin flip-up wheels



Convertible armrests

FEATURES

- Convertible Armrests are both adjustable in height and length from Desk to Full in three settings
- Stylish rear Mag Wheels have brushed aluminum handrims and quick release rear axles for fast, too free hemi-height adjustment and for easy breakdown and transport
- Comes included with a seat rail and upholstery extension that easily adjusts seat depth from 16" to 18" without any additional upholstery
- The chair comes standard with a wrap-around upholstery and an adjustable height and angle back (from 5 degrees to 20 degrees)
- Comfortable, slick nylon back upholstery is fluid resistant and easy-to-clean
- Adjustable anti-tippers included
- HCPCS* Code: K0004

K0195 or E0990 (elevating legrests only)
 E0971 (anti-tippers)
 E2201 (non-standard seat frame - 20" chair only)

SPECIFICATIONS

| | |
|-------------------------------------|--|
| Seat Width: | 16" 18" 20" 22" |
| Seat Depth: | 16" - 18" |
| Seat-to-Floor Height ¹ : | Adult 19.5" Hemi 17.5" |
| Back Height: | 16" - 18" |
| Back Style: | Height Adjustable |
| Overall Width w/Arms: | 24.5" 26.5" 28.5" |
| Overall Depth: | 32" (without riggings) |
| Weight Capacity: | 300 LBS (16" / 18" / 20") 350 LBS (22") |
| Product Weight ² : | 32 - 36 lb |
| Carton Dimensions: | 33" x 13" x 37" |
| Shipping Weight: | 35 - 39 lb |
| Warranty: | Limited Lifetime on Frame, One Year on all Wear Parts |

WARNING: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer or birth defects or other reproductive harm. For more information, visit www.P65Warnings.ca.gov

| Item # | Seat | Arm Style | Leg Style | UOM |
|-------------|-----------|---|----------------------|------|
| WC4DX1616CS | 16" x 16" | Height Adjustable, Flip Back, Convertible | Swing-Away Footrests | each |
| WC4DX1816CS | 18" x 16" | Height Adjustable, Flip Back, Convertible | Swing-Away Footrests | each |
| WC4DX2016CS | 20" x 16" | Height Adjustable, Flip Back, Convertible | Swing-Away Footrests | each |
| WC4DX2216CS | 22" x 16" | Height Adjustable, Flip Back, Convertible | Swing-Away Footrests | each |
| WC4DX1616CE | 16" x 16" | Height Adjustable, Flip Back, Convertible | Elevating Legrests | each |
| WC4DX1816CE | 18" x 16" | Height Adjustable, Flip Back, Convertible | Elevating Legrests | each |
| WC4DX2016CE | 20" x 16" | Height Adjustable, Flip Back, Convertible | Elevating Legrests | each |
| WC4DX2216CE | 22" x 16" | Height Adjustable, Flip Back, Convertible | Elevating Legrests | each |

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary 1/4-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.