



K4 Deluxe Wheelchair

HIGH STRENGTH K4 WHEELCHAIR

The **ProBasics K4 Deluxe Wheelchair** is an high-strength, feature rich wheelchair. This wheelchair comes standard with a 300 LB weight capacity, One -Button-Convertible Desk-To-Full Length Armrests, Quick Release Axle for Hemi-Height Adjustments. The K4 Deluxe easily adapts to any patient's individual needs. It offers superior value and functionality for those looking for a high-performance mobility product.



2" upholstery seat extension



Adjustable anti-tippers with push-pin flip-up wheels



Quick release rear axle

FEATURES

- Single Button Armrest Adjustment. Single Button to change from desk length to full length arms.
- Stylish rear Mag Wheels have brushed aluminum handrims and quick release rear axles for fast, too free hemi-height adjustment and for easy breakdown and transport.
- Comes included with a seat rail and upholstery extension that easily adjusts seat depth from 16" to 18" without any additional upholstery
- Adjustable height/angle wrap-around back upholstery for maximum comfort.
- Comfortable, slick nylon back upholstery is fluid resistant and easy-to-clean.
- Adjustable anti-tippers included
- HCPCS* Code: K0004

K0195 or E0990 (elevating legrests only)

E0971 (anti-tippers)

E2201 (non-standard seat frame - 20" chair only)

SPECIFICATIONS

Seat Width:	16" 18" 20" 22"
Seat Depth:	16" - 18"
Seat-to-Floor Height ¹ :	Adult 19.5" Hemi 17.5"
Back Height:	16" - 18"
Back Style:	Height Adjustable
Overall Width w/Arms:	24.5" 26.5" 28.5"
Overall Depth:	32" (without riggings)
Weight Capacity:	300 lb (16/18/20") / 350 lb. (22")
Product Weight ² :	32 - 36 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	35 - 39 lb
Warranty:	Limited Lifetime on Frame, One Year on all Wear Parts

Item #	Seat	Arm Style	Leg Style	UOM
WC4DX1616CS	16" x 16"	Height Adjustable, Flip Back, Convertible	Swing-Away Footrests	each
WC4DX1816CS	18" x 16"	Height Adjustable, Flip Back, Convertible	Swing-Away Footrests	each
WC4DX2016CS	20" x 16"	Height Adjustable, Flip Back, Convertible	Swing-Away Footrests	each
WC4DX2216CS	22" x 16"	Height Adjustable, Flip Back, Convertible	Swing-Away Footrests	each
WC4DX1616CE	16" x 16"	Height Adjustable, Flip Back, Convertible	Elevating Legrests	each
WC4DX1816CE	18" x 16"	Height Adjustable, Flip Back, Convertible	Elevating Legrests	each
WC4DX2016CE	20" x 16"	Height Adjustable, Flip Back, Convertible	Elevating Legrests	each
WC4DX2216CE	22" x 16"	Height Adjustable, Flip Back, Convertible	Elevating Legrests	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary 1/4-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.