



## ProBasics Reclining Manual Wheelchair

The **ProBasics Reclining Manual Wheelchair** provides full back and headrest support as patient reclines from 90 to 160 degrees. The 13-inch headrest and padded, full-length arms are removable. Includes black, padded vinyl upholstery. The frame has a chip-resistant silver vein finish. Mag-style wheels are lightweight and set back on frame to prevent patient from tipping when in the reclined position.



Reclines from 90° – 160°



Caregiver recliner adjustment mechanism



13-inch headrest supports patient when reclined and is removable

### FEATURES

- Slide tube steel frame with an attractive silver vein finish that is chip resistant and easy to maintain
- Padded Desk-length, detachable armrests provide added patient comfort
- Easy-to-clean, comfortable vinyl upholstery
- Chart pocket on backrest provides additional convenience
- Tool-free elevating padded legrests are equipped with a pneumatic piston for infinite positions and composite footplates
- Composite handrims with mag-style wheels are lightweight and maintenance free
- Wheels are set back on frame to prevent tipping
- Rear anti-tippers come standard
- Headrest extension comes standard with head immobilizer
- 8" x 1" solid rubber front casters
- Easy push-to-lock wheel locks
- HCPCS\* Code: K0001/K0006

E1226 (reclining back)

E0966 (headrest extension)

E0990 (elevating legrests only)

### SPECIFICATIONS

Seat Width:	16"   18"   20"   22"
Seat Depth:	17"
Seat-to-Floor Height <sup>1</sup> :	Adult 20"
Back Height:	19.5"   34"
Back Style:	Reclining
Overall Width w/Arms:	24½"   26½"   28½"
Overall Depth:	31¼" (without riggings)
Weight Capacity:	300 lb   400 lb (WCR2218E)
Product Weight <sup>2</sup> :	51 – 60 lb
Carton Dimensions:	40" x 12" x 36"
Shipping Weight:	71 – 75 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	Weight Limit	UPC	UOM
WCR1616E	16" x 17"	Removable, Desk-Length	Elevating Legrests	300 lb	815067070893	Each
WCR1816E	18" x 17"	Removable, Desk-Length	Elevating Legrests	300 lb	815067070909	Each
WCR2018E	20" x 17"	Removable, Desk-Length	Elevating Legrests	300 lb	815067070916	Each
WCR2218E	22" x 17"	Removable, Desk-Length	Elevating Legrests	400 lb	815067070923	Each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

<sup>1</sup>NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

<sup>2</sup>NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.