



ProBasics K4 Transformer

TRANSFORMS FROM K4 TO A TRANSPORT CHAIR

The perfect upgrade from a standard K0004 wheelchair, the **ProBasics TRANSFORMER™** is a dual purpose, high strength lightweight wheelchair that quickly transforms from a self-propelled wheelchair to a compact transport chair with 8" wheels all around. Great as a transitional chair for accompanied excursions while retaining self-propel capabilities. Also allows the slimming down of the chair to navigate narrow doorways.



24-inch quick release rear wheels



Transforms into a Transport Chair



Flip back, desk length, padded arms

FEATURES

- Lightweight aluminum frame in blue flame color
- 8-inch full swivel front casters mounted on aluminum forks
- 24-inch "quick release" rear wheels with composite handrims
- Bolted on 8-inch transport chair rear wheels
- Self-locking paddle release, flip-back padded desk length arms
- Lightweight black, woven nylon upholstery
- Fold down back with easy release paddles
- Complete with lightweight aluminum wheel locks for both 24-inch and 8-inch rear wheel modes
- Detachable flip-up rear anti-tippers
- Accessory pouch sewn into back
- Footrests include heel loops and composite footplates (optional elevating legrests sold separately - WCT4ELR)
- Folds for easy storage
- HCPCS* Code: K0004

K0195 or E0990 (elevating legrests only)

SPECIFICATIONS

Seat Width:	16" 18" 20"
Seat Depth:	16"
Seat-to-Floor Height ¹ :	Adult 20"
Back Height:	16"
Back Style:	Folds downs
Overall Width w/Arms:	18" - 29"
Overall Depth:	28" - 32" (without riggings)
Weight Capacity:	250 lb
Product Weight ² :	32 - 34 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	35 - 37 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WCT41616DS	16" x 16"	Flip Back, Desk-Length	Swing-Away Footrests	815067070817	each
WCT41816DS	18" x 16"	Flip Back, Desk-Length	Swing-Away Footrests	815067070824	each
WCT42016DS	20" x 16"	Flip Back, Desk-Length	Swing-Away Footrests	815067070831	each
Item #	Description			UPC	UOM
WCT4ELR	Optional Elevating Leg Rests			815067070848	pair

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.