


Patient Interfaces

Sleepnet® Aura™ CPAP Mask


Enjoy comfort and a great seal with the world's first nasal mask featuring both CustomFit Technology and an Active Headgear Connector. With three sizes and a Flexible Spacebar, Aura will give your patients the rest they deserve.

 This product is NOT manufactured with natural rubber latex.



Product Specifications

- Soft shell mask has a shapeable wire allowing for optimum flexibility, fit and comfort
- Active Headgear Connector moves when patient moves to help maintain a seal throughout the night
- AIR°gel with Advanced Cushion Technology creates an ergonomic gel cushion that reduces pressure points on the skin
- Available in three sizes with simple and breathable quick release EZ Fit headgear
- Made in the USA



Item #	Description	Qty
50176	Sleepnet® Aura™ CPAP Mask, Small	each
50177	Sleepnet® Aura™ CPAP Mask, Medium	each
50178	Sleepnet® Aura™ CPAP Mask, Large	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com