

# Cliq™ Aspirator

## VC-701 Model

A single-cylinder aspirator producing high-vacuum, low-flow for medical suction. Portable and durable design. Includes 800 cc collection canister with overflow protection, microbial filter with ¼" hose barb, 2 silicone intermediate tubes and 6' suction tubing. Suction cups leg mounts provide extra stability while in use.

### Product Specifications

Max Air Flow	-15-17 LPM
Max Vacuum Pressure	-675 mmHg
Pressure Display	Analog Gauge
Weight	7.7 lbs
Dimensions	13" x 6" x 8.25"
Electrical Requirements	110-120V/220-240V (50-60Hz)
Power Cord Length	6'
Operating Temperature	32° - 104° F
Operating Relative Humidity	0-95%
Storage and Transport Temp	-40° - 158° F
Storage and Transport Humidity	0-95%
Warranty	Two Years
HCPCS Code	E0600



Item #	Description	Qty	Case QTY
CF608	Cliq™ Aspirator	Each	1
ILB-ROS	Suction Filter	Each	1
50038	800 cc Collection Canister with Floater Top	Case	12
TK-MA	Silicone Tubing, 7"	Pack	10