

# Cliq™ 50 psi Compressor

## CX-500 Model

The CX-500 is a 50 psi compressor designed for efficient high pressure performance for extended periods of time. It is intended to administer humidified air or medication into the airways to treat respiratory disorders. Designed for use with heated nebulizers, compatible hand-held nebulizer systems, dual action humidifiers/nebulizers, mist tents, high output pneumatic nebulizers and deluxe nebulizers with oxygen diluter. Includes U-shaped DISS connection for single neb bottle and 2 extra air filters.

### Product Specifications

Air Flow	0 psi: 37 LPM 50 psi: 17 LPM
Max Pressure	95 psi (preset)
Pressure Display	Analog Gauge
Weight	14.9 lbs
Dimensions	14" x 7.5" x 12.5"
Electrical Requirements	110-120V / 220-240V (50-60 Hz)
Operating Temperature Range	41° - 104° F
Operating Relative Humidity	0-93%
Storage and Transport Temperature	-13° - 158° F
Storage and Transport Humidity	0-93%
Warranty	One Year
HCPCS Code	E0565



Item #	Description	Qty	Case QTY
CX-500	Cliq™ 50 psi Compressor	Each	1
HCS00350	Prefilled Humidifier, 350 mL, Sterile Water for Inhalation USP	Case	20