

## Medical Gas

# Adjustable Oxygen Carts

Roscoe Adjustable Oxygen Carts are constructed of durable, high-strength steel and can accommodate either "D" or "E" size cylinders. Handle is height adjustable to accommodate most individuals. Large, easy-roll wheels allow for safe maneuvering. Includes gurney/mask hook, nonslip hand grip and one piece thumb screw. Available in silver vein finish or with chrome plating. Limited lifetime warranty. HCPCS Code: E1355\*



E-CART



E-CARTPC



E2-CARTPC

Item #	Description	Height Adjustable	Cylinder Capacity	Weight	Qty
E-CART	D or E Cylinder Cart, Chrome Finish	38 1/2" - 43"	1	5 lb.	4/cs
E-CARTPC	D or E Cylinder Cart, Silver Vein Finish	38 1/2" - 42 1/2"	1	5 lb.	4/cs
E2-CARTPC	D or E Cylinder Cart, Silver Vein Finish	38 1/2" - 42"	2	9 lb.	2/cs

## Replacement Parts for Oxygen Cart



E-SCREW10



90181



90183

Item #	Description	Qty
E-FOOTPK	Rubber Feet for E-CART, E-CARTPC, E2-CARTPC	2/pk
E-SCREW10	Thumb Screw for E-CART	10/pk
90181	Rubber Grip, Black for E-CART and E-CARTPC	each
90185	Rubber Grip, Black for E2-CARTPC	each
90183	Replacement Wheel for E-CART	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



**To Order:**  
 (800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com