

## Electrotherapy

### EMS 5.0 Muscle Stim

The EMS 5.0 Muscle Stim Analog Unit is designed for personal use, however it is strong enough to be used in an office setting. The EMS unit uses electrical impulses to contract and relax muscles helping those who experience muscle tension in the upper back and neck areas. Can be used to regain muscle strength after surgery or being in a cast. The EMS unit is portable and may be clipped to a belt, shirt, pocket or other clothing.



#### Product Features:

- System includes device, lead wires, four self-adhesive reusable electrodes, 9-volt battery, hard plastic carrying case and instruction manual
- 105 mA strong
- Power Source: 9-volt battery
- Channel: Dual, isolated between channels
- Pulse Width: 250 microseconds
- Pulse Rate: 5, 30, 100 Hz
- Pulse Ramp: 1, 3, 5 seconds
- Pulse Amplitude: Adjustable, 0-80 mA
- Contraction Control: 1-30 seconds
- Relaxation Control: 1-45 seconds
- Wave Form: Asymmetrical Bi-Phase Square Pulse
- One year warranty

Item #	Description	Qty
DE5030	EMS 5.0 Muscle Stim Analog Unit	each

### EMS 7500 Digital Edition

The EMS 7500 Digital Device features a patient compliance meter, three modes, safety amplitude cover, sturdy metal belt clip, rubber side railings for maximum handling and a premium silver finish.



#### Product Features:

- System includes device, lead wires, four self-adhesive reusable electrodes, 9-volt battery, hard plastic carrying case and instruction manual
- Patient compliance meter can store up to 60 sets of operation records and 999 hours
- Power Source: 9-volt battery
- Channel: Dual, isolated between channels
- Pulse Width: 50-300  $\mu$ s, adjustable
- Pulse Rate: 2-120 Hz, adjustable
- Pulse Ramp: 1-8 seconds
- Pulse Amplitude: 100 mA peak into 500 phm load
- Contraction Time: 0-90 seconds
- Relaxation Time: 2-90 seconds
- Waveform: Asymmetrical Bi-Phase Square Pulse
- Timer: 1-60 minutes or continuous
- One year warranty

Item #	Description	Qty
DE7502	EMS 7500 Digital Edition	each

\* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



**To Order:**

(800) 3-ROSCOE  
fax: (440) 572-4261  
roscoemedical.com

## Electrotherapy

# Twin Stim® TENS and EMS Combo

Our most popular digital electrotherapy unit. The Twin Stim 2nd Edition is an exact combination of our TENS 7000 and EMS 7500, which gives the user the choice of five Modes of TENS and two Modes of EMS. Features include, dual channels, seven modes, timer, patient compliance meter, new stronger 100 mA output, safety amplitude cover, sturdy metal belt clip and rubber side railings for maximum handling.



### Features:

- Five TENS Modes: Burst / Normal / Modulation / SD1 / SD2
- Two EMS Modes: Synchronous, Alternate
- System includes device, lead wires, four self-adhesive reusable electrodes, 9-volt battery, hard plastic carrying case and instruction manual
- Patient compliance meter can store up to 60 sets of operation manual
- Power Source: 9-volt battery
- Pulse Width: 50-300  $\mu$ s, adjustable
- Pulse Rate: 2-150 Hz, adjustable
- Pulse Ramp: 1-8 seconds
- Pulse Amplitude: 100 mA peak into 500 ohm load
- Contraction Time: 2-90 seconds
- Relaxation Time: 2-90 seconds
- Wave Form: Asymmetrical Bi-Phase Square Pulse
- Timer: 1-60 minutes or continuous
- One year warranty

Item #	Description	Qty
DS2202	Twin Stim® TENS and EMS Combo 2nd Edition	each

## Twin Stim® Plus

The popular Twin Stim® Plus is a four channel TENS and EMS combo unit complete with eight electrodes. The Twin Stim® Plus gives the user the choice of five modes of TENS and three modes of EMS. Unit can be powered by four AA batteries or by AC adaptor.



### Product Features:

- Five TENS Modes: Burst / Normal / Modulation / SD1 / SD2
- Three EMS Modes: Constant / Synchronous, / Alternate /
- System includes device, AC adaptor, four lead wires, eight self-adhesive reusable electrodes, four AA batteries, hard plastic carrying case and instruction manual
- Built in patient compliance monitor
- Power Source: 4 AA batteries
- Pulse Width: 50-300  $\mu$ s, adjustable
- Pulse Rate: 2-150 Hz, adjustable
- Pulse Ramp: 1-8 seconds
- Pulse Amplitude: Adjustable: 100 mA peak into 500 ohm load
- Contraction Time: 2-90 seconds, 1 sec./step
- Relaxation Time: 0-90 seconds, 1 sec./step
- Wave Form: Asymmetrical Bi-Phase Square Pulse
- Timer: 1-60 minutes or continuous
- One year warranty

Item #	Description	Qty
DS5002	Twin Stim® Plus, Four Channel TENS and EMS Combo	each

\* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



**To Order:**

(800) 3-ROSCOE  
fax: (440) 572-4261  
roscoemedical.com