

# PORTABLE ELECTROTHERAPY DEVICES FOR PAIN RELIEF

InTENSity portable electrotherapy devices provide the most up-to-date technology and waveforms (TENS, NMES, IF, and/or Russian), with preset parameters for use on specific areas of the body. Each device provides 2 channels of output, and is powered by a rechargeable lithium ion battery to be used at home, work, or on the go.



## INTENSITY™ 12 REF DI1212

- TENS waveform
  - 4 TENS modes:
    - continuous
    - pulse width modulation
    - pulse rate modulation
    - combo modulation
- 12 Preset Modes & 6 Manual User Settings
- Preset and adjustable treatment times
- Store up to 90 sets of treatment records



## INTENSITY™ TWIN STIM IV REF DI2717

- TENS and NMES
  - 4 TENS modes:
    - continuous
    - pulse width modulation
    - pulse rate modulation
    - combo modulation
  - Synchronous mode for NMES with adjustable parameters:
    - pulse amplitude
    - pulse frequency
    - contraction/relaxation time
- 24 Preset Modes & 12 Manual User Settings
- Preset and adjustable treatment times
- Store up to 90 sets of treatment records



## INTENSITY™ IF COMBO II REF DI2738

- TENS and Interferential Stimulation
  - 4 TENS modes:
    - continuous
    - pulse width modulation
    - pulse rate modulation
    - combo modulation
  - Up to 10,000 Hz of sweeping IF to prevent accommodation of nerves and provide maximum treatment benefits
- 24 Preset Modes & 12 Manual User Settings
- Preset and adjustable treatment times
- Store up to 90 sets of treatment records



## INTENSITY™ SELECT COMBO II REF DI2195

- TENS, NMES, IF, and Russian waveforms
  - 4 TENS modes:
    - continuous
    - pulse width modulation
    - pulse rate modulation
    - combo modulation
  - Synchronous mode for NMES with adjustable parameters:
    - pulse amplitude
    - pulse frequency
    - contraction/relaxation time
  - Up to 10,000 Hz of sweeping IF to prevent accommodation of nerves and provide maximum treatment benefits
  - 2500 Hz of Russian waveform for deep and comfortable stimulation
- 48 Preset Modes & 24 Manual User Settings
- Preset and adjustable treatment times
- Store up to 90 sets of treatment records



Each InTENSity kit includes the portable stimulation device, a pack of four electrodes, a pair of lead wires, 1 lithium ion battery, a USB wall charger, belt clip, instruction manual, quick start guide, and carrying case. InTENSity portable electrotherapy devices are covered by a one year warranty.