



Roscoe Knee Scooter

The **Roscoe Knee Scooter** offers increased comfort and ease of mobility to patients who cannot put weight on their foot or ankle, but want to remain active. With its four large wheels and padded knee platform, the Knee Scooter is perfect for indoor and outdoor use. The offset post on the aluminum knee platform allows for right-to-left adjustment in seconds - no tools required. The Knee Scooter is equipped with an easy folding mechanism that collapses using a thumb release lever and fits beneath the raised knee platform for ease of storage and transportation.



Folds Flat for Easy Storage



Adjustable Hand Brake

FEATURES

- Scooter can be steered for increased maneuverability
- Ideal for individuals recovering from an ankle or foot injury
- Large 8-inch clear PU wheels with sealed bearings
- Adjustable hand brake ensures safety
- Padded knee platform with easy-adjustment mechanism
- Padded handlebar grips
- Limited lifetime warranty on frame, 90 days on all other parts

COLORS

Black
Finish



White
Finish



Burgundy
Finish



SPECIFICATIONS

Patient Height Range: 4'11" – 6'6"

Knee Platform Dimensions: 7" (W) x 14" (L)

Adjustable Handle Height: 31" – 41"

Knee Platform Height: 17.5" – 22"

Product Weight: 22.8 lbs.

Weight Capacity: 350 lbs.

Item #	Description	UPC	UOM
ROS-KSB	Roscoe Black Knee Scooter	092237601136	each
ROS-KSW	Roscoe White Knee Scooter	092237601471	each
ROS-KSBG	Roscoe Burgundy Knee Scooter	092237601143	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.