

Roscoe Mobility

Kona KI/K2 Wheelchair

The Roscoe Kona™ Dual Axle KI/K2 Manual Wheelchair combines style and durability at an affordable price. Standard features include a dual axle, comfortable nylon upholstery, padded armrests, chart pocket and front riggings. The Kona™ Wheelchair supports patient weights up to 250 pounds with the option of swing-away footrests or elevating legrests.



Silver Vein Finish



Product Features

- Slide tube, powder-coated silver vein steel frame
- Dual axle positioning
- Easy-to-clean, comfortable nylon upholstery resists mildew and bacteria
- Padded and upholstered armrests come standard
- Chart pocket on backrest provides additional convenience
- Includes tool-free swing-away footrests or elevating legrests with composite footplates
- Composite handrims with 24" solid rubber tires
- 8" x 1" solid rubber front casters
- Easy push-to-lock wheel locks
- 5 year warranty on frame and cross braces
90 days on all other parts

Product Specifications

- Product Weight²: 38 lb.
- Overall Width with Arms: 26 ½"
- Overall Depth: 31" (without front riggings)
- Seat Width: 16" / 18" / 20"
- Seat Depth: 16"
- Back Height: 18"
- Seat-to-Floor Height¹: Adult 20" / Hemi 18"
- Back Style: Fixed
- Weight Capacity: 250 lb.
- HCPCS Code: K0001, K0002*



Composite Footplates



Chart Pocket



Push-Style Wheel Locks



Wheelchair Replacement Parts Available!

Item #	Frame	Seat Width	Armrests	Front Rigging
W11616DS	Slide Tube	16"	Removable Desk-Length Arms	Swing-Away Footrests
W11616DE	Slide Tube	16"	Removable Desk-Length Arms	Elevating Legrests
W11816DS	Slide Tube	18"	Removable Desk-Length Arms	Swing-Away Footrests
W11816DE	Slide Tube	18"	Removable Desk-Length Arms	Elevating Legrests
W12016DS	Slide Tube	20"	Removable Desk-Length Arms	Swing-Away Footrests
W12016DE	Slide Tube	20"	Removable Desk-Length Arms	Elevating Legrests

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.
¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.
²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how wheelchair is equipped.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com