

Patient Interfaces

Sleepnet® Mojo® Full Face Mask

Sleepnet's® Mojo® Full Face Mask features a contoured shell that follows the curves of the human face. Air leakage is eliminated by adjusting the headgear, as the dual beam element transfers the force of the straps to the center of the mask. Sleepnet's® unique magnetic clips are virtually foolproof, and easily click into position securing the mask and keeping the straps in place. The clips allow for a quick release of the mask and are designed to withstand the rigors of lab or home use.



This product is NOT manufactured with natural rubber latex.

Product Specifications

- Soft-shell mask has a shapeable wire allowing for optimal flexibility, fit and comfort
- The Adjustable Tension Wheel is exclusive to the Sleepnet® Mojo® - a wide, continuous range of adjustment ensures a custom fit and removes pressure from the bridge of the nose.
- Available in three sizes with magnetic connectors and simple, breathable quick-release headgear

AIR°gel

Maximum comfort, minimum weight. Our patented AIR°gel minimizes pressure points and skin irritation.



AIR°gel Forehead Cushion

Adjustable Tension Wheel

Removes pressure from the bridge of the nose for a custom, comfortable fit.

Custom Fit Technology™

The shapeable shell of each mask can be easily molded to fit individual needs. **Flex, fit, forget.**

Easy-Release Magnetic Headgear Clips

Lightweight, breathable headgear disconnects from mask easily with one hand.

360° Rotating Elbow

For optimal tube positioning and easier mobility.

Vented and Non-Vented Versions Available



Item #	Description	Qty
50835	Sleepnet® Mojo® Full Face Mask, Small	each
50836	Sleepnet® Mojo® Full Face Mask, Medium	each
50837	Sleepnet® Mojo® Full Face Mask, Large	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:

(800) 3-ROSCOE
fax: (440) 572-4261
roscoemedical.com