

## Home Care Beds

### Non-Tilt Overbed Table

The Roscoe Overbed Table is practical, easy-to-use and suitable for both home and institutional use. It is versatile enough for eating, reading and many other patient activities. Our Overbed Table is easy to assemble and comes equipped with swivel casters to provide smooth rolling action. The attractive walnut wood-grain laminate tabletop is designed to resist everyday wear-and-tear.



#### Product Features

- Can be easily raised or lowered in infinite positions
- Chrome-plated steel "H" base provides security and stability
- Swivel casters allow for easy maneuverability
- Overall width at base is 26 ½ inches
- Table top can be raised with the slightest upward pressure
- Weight capacity: 40 lb.
- Limited lifetime warranty
- HCPCS Code: E0274\*

Overbed Table Accessories	Item #	Qty
Replacement Table Top	90011	1
Replacement Casters	90009	4
Hardware Set (blister pack)	90007	1

Item #	Tabletop Adjustment (highest)	Tabletop Adjustment (lowest)	Tabletop Surface (width)	Tabletop Surface (depth)	Product Weight	Qty
ROS-OBT	47"	30"	30"	15"	22 lb.	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



**To Order:**  
 (800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com