

Tracheal Care

Roscoe Portable Suction Machine

The lightweight Roscoe Medical Portable Suction Machine is easy to handle and transport, perfect for everyday use. The compact design has a fully enclosed handle that offers patients a less institutional look. The unit is available in AC power or with DC rechargeable battery power. Designed for easy portability at just 5.7 pounds, the Roscoe Portable Suction Machine includes a 1000 cc collection bottle and a two year warranty.



Product Features

- Available in AC power or with DC rechargeable battery
- Compact design, perfect for everyday use
- Handle is fully integrated offering a more contemporary look
- Complete kit includes aspirator, reusable collection jar, filter with 6' suction tubing, carry bag and instruction manual
- Rechargeable lithium battery charges in less than eight hours (50006 only)
- The working cycle of the device is 30 minutes on and 30 minutes off
- 90 day warranty on battery
- Two year warranty on device
- HCPCS Code: E0600*

Technical Specifications

- Dimensions: 14.3" (L) x 6.9" (W) x 8.3" (H)
- Weight: 5.7 lb. (50004), 5.3 lb. (50006)
- Maximum Vacuum: 550 mmHg (50004), 530 mmHg (50006)
- Flow Range: 28 LPM (50004), 25 LPM (50006)
- Sound Level: ≤ 65 dBA
- Electrical Requirements: 120V/60Hz (50004), 12V (50006)
- Power Consumption: 70 watt (50004), 48 VA (50006)
- Collection Capacity: 1000 cc
- Operating Temperature: 50° - 104°F
- Operating Humidity: 95% RH

WARNING: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer or birth defects or other reproductive harm. For more information, visit www.P65Warnings.ca.gov

Item #	Description	Qty
50004	Roscoe Medical Portable Suction Machine with AC Power	each
50006	Roscoe Medical Portable Suction Machine with DC Rechargeable Battery	each

Replacement Parts and Accessories

Item #	Description	Qty
50038	800 cc Disposable Collection Jar	12/pk
KIT-50004	Suction Kit (includes elbow, bacteria filter, collection jar, 6' tubing and two 7" tubing)	1/kt
90107	Battery Charger (Brick)	each
90109	Gauge	each
90111	Power Cord	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
roscoemedical.com