

## Alternating Pressure Systems

### Quick-Fill Pump/Magic-Air®

The Quick-Fill Pump is designed to quickly inflate and deflate Low Air Loss Mattresses reducing your patient set-up time by 65%, and allowing you time to attend to your patient's other needs. This quick 20-liter-per-minute flow pump allows your technician to get the mattress up and therapeutic in minutes.



### Product Features

- Light weight, compact and all-steel construction
- Pumps air at a rate of 20-liters-per-minute
- Quieter as compared to standard plastic air pumps
- Weight: 3 lbs.
- Dimensions: 13" length x 5.5" wide height
- 10 foot cord
- Works with Roscoe Medical's Quick Connect Adaptor 90121
- One year warranty

\*Magic-Air® is a registered trademark of Metropolitan Vacuum Cleaner Co., Inc.

Item #	Description	Qty
90501	Quick-Fill Pump/Magic-Air®	each
90153	Quick Fill Valve	each



**To Order:**  
 (800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com