


Medical Gas

Cylinder Racks and Stands

Roscoe Medical's Cylinder Racks and Stands feature an all-welded frame construction to ensure durability and strength. Individual spacing, positive cylinder placement and weight distribution promote safe storage. Equipped with bolt down bars and feet, our stands are available in chrome or a polyurethane rubber-coated green finish, designed to protect and preserve the surface of your aluminum cylinders.




CR-M6PC6



CR-M6R6



CR-M6PC12



CR-M6R12



CR-MST

Item	Description	Finish	Capacity	Size (inches)	Weight
CR-M6PC6	M6 Cylinder Rack, 6 hole	Chrome	6 Cylinders	12.8 x 8.46 x 9.45	5 lbs.
CR-M6R6	M6 Cylinder Rack, 6 hole	Green, PU Coated	6 Cylinders	12.8 x 8.46 x 9.45	5 lbs.
CR-EPC6	C/D/E Cylinder Rack, 6 hole	Chrome	6 Cylinders	16.34 x 10.83 x 9.65	8 lbs.
CR-ER6	C/D/E Cylinder Rack, 6 hole	Green, PU Coated	6 Cylinders	16.34 x 10.83 x 9.65	8 lbs.
CR-M6PC12	M6 Cylinder Rack, 12 hole with crossbars	Chrome	12 Cylinders	17.72 x 12.8 x 9.45	9 lbs.
CR-M6R12	M6 Cylinder Rack, 12 hole with crossbars	Green, PU Coated	12 Cylinders	17.72 x 12.8 x 9.45	9 lbs.
CR-EPC12	C/D/E Cylinder Rack, 12 hole with crossbars	Chrome	12 Cylinders	21.65 x 17.72 x 10.24	14 lbs.
CR-ER12	C/D/E Cylinder Rack, 12 hole with crossbars	Green, PU Coated	12 Cylinders	21.65 x 17.72 x 10.24	14 lbs.
CR-MST	M60, M90 and MM Cylinder Stand	Green, Powder Coated	1 Cylinder	14.57 x 13.78 x 11.22	5 lbs.



To Order:

(800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com