

Aerosol Therapy

Roscoe SierraNeb 2 Nebulizer System

The Roscoe SierraNeb 2 Nebulizer System offers superior performance at an affordable price. The simplistic one-button operation and ultra-compact lightweight design makes it extremely portable and patient friendly. Weighing only 2.65 pounds, the SierraNeb 2 is easy to handle and transport, providing reliable treatment anywhere. The complete system includes compressor, nebulizer, angled mouthpiece, tubing, five extra filters, instruction manual and five year warranty. The Roscoe SierraNeb 2 Nebulizer System is also available with a reusable neb kit option.

Product Specifications

- One-button operation
- Powerful piston pump is efficient, dependable and quiet
- Ultra compact design makes it extremely portable
- Produces a particle size of less than 5 microns
- Darker body color maintains its appearance after prolonged periods of use
- Comes complete with compressor, nebulizer, angled mouthpiece, tubing, five extra filters and instruction manual
- TruNeb Reusable Neb Kit available - purchase with the compressor or alone
- Attractive retail packaging
- Five year warranty
- HCPCS Code: E0570*
- HCPCS Code: A7003* for nebulizer kits
- HCPCS Code: A7005* for reusable nebulizer kits
- HCPCS Code: A7013* for foam filters



NEB-SIERRA

Technical Specifications

- Power: 120 VAC/60Hz
- Power Consumption: 1.4 amps
- Sound Level: 58 dBA at 1 meter
- Maximum Compressor Pressure: 35 - 45 PSI
- Typical Operating Flow Rate: 12 PSI
- Average Nebulization Rate: 0.25 - 0.35 ml/min
- Medication Capacity: 5 ml
- Particle Size: Less than 5 microns
- Operating Temperature: 50° - 104° F
- Operation Cycles: 20 minutes on, 40 minutes off
- Dimensions: 6 (L) x 4" (W) x 3" (H)
- Weight: 2.65 lb.

Accessories
Included

BAG-NEBR



BAG-DRW

Item #	Description	Qty
NEB-SIERRA	Roscoe SierraNeb 2 Nebulizer System with Disposable Neb Kit	each
SIERRA-TRU	Roscoe SierraNeb 2 Nebulizer System with one Disposable Neb Kit and one TruNeb Reusable Neb Kit	each
BAG-NEBR	Nylon Carry Bag for SierraNeb 2 Nebulizer System	each
BAG-DRW	Nylon Drawstring Carry Bag for Mini Nebulizer Compressor System, 13" x 12"	each
CIF-1128PK6	Foam Filter, 5/8" diameter x 1/2" thick	6/pk



* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

To Order:

(800) 3-ROSCOE
fax: (440) 572-4261
roscoemedical.com