

truAIRE-PASP Portable Aspirator

The truAIRE-PASP Portable Aspirator is engineered for rapid response and reliable performance, the truAIRE-PASP delivers powerful suction in a compact, lightweight design. Ideal for emergency medical services and home healthcare, this unit offers exceptional portability without compromising suction strength. With a simple flow control adjustment knob and long-lasting rechargeable battery, this device is ready to go when it matters the most.



Easy to read vacuum gauge and variable suction settings



True portability - Unit can operate while inside carry bag



Unit can be charged while in use

PRODUCT FEATURES:

- Compact and portable design, perfect for oral, nasal, and tracheal suctioning
- Can be operated while in carry bag for on-the-go use
- Variable vacuum regulation from 50 to 550 mmHg
- Complete with carry bag, dual elbow tubing assembly, 800cc collection jar, bacteria filter, 6' suction tubing, and user manual
- Replaceable, lithium battery has a 90-minute duration with a working cycle of 30 minutes on and 30 minutes off
- Comes complete with A/C D/C power supply with 90° connector
- 2 Year warranty on the device, 1 Year warranty on the battery
- HCPCS Code: E0600*

TECHNICAL SPECIFICATIONS:

Dimensions:	8.25" L x 7.5" W x 7.75" H
Weight:	3.6 lbs
Maximum Vacuum:	600mmHg (Recommend Range 50 – 550mmHg)
Flow Range:	18 LPM +/- 10%
Sound Level:	≤ 60 dBA
Electrical Requirements:	100-240V AC, 50-60Hz
Operating Temperature:	32° - 104°F (0°-70°C)
Operating Humidity:	0 - 95%

WARNING: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer or birth defects or other reproductive harm. For more information, visit www.P65Warnings.ca.gov

Item #	Description	UPC	Qty
ASP-ROSP	truAIRE-PASP Portable Aspirator	092237626672	2/CA
ASP-KIT	Disposable Suction Kit for ASP-ROSP (Patient Kit contains Canister, 2-Elbow Tubes, 6' Tube & Bacteria Filter)		1/ea
1160B-1	Bemis 800cc Suction Canister & LID		70
55N56A	6' Suction Tubing		50
X1PASP01	Elbow Tubing Assembly for ASP-ROSP		1/ea
X1PASP03	Bacteria Filter for ASP-ROSP		1/ea
PORT-CARADP	DC Car Charger		1/ea
X1PASP04	Replacement Battery - ASP-ROSP		1/ea