



CPAP CATALOG





CPAP PRODUCTS

Roscoe® has consistently offered tools that help sleep professionals provide better care and achieve better results. We are committed to help you break through emerging clinical and industry barriers by offering more options in CPAP Therapy than ever before. Find everything you need at Roscoe to manage your patients and your business.

CPAP Products

| | |
|---------------------------|-------|
| Patient Interfaces | 4–8 |
| CPAP Systems | 9 |
| CPAP Private Label | 10 |
| CPAP Filters | 11–12 |
| CPAP Tubing | 13 |
| CPAP Replacement Headgear | 15 |
| CPAP Accessories | 16–21 |

DreamEasy™ 2 Full Face CPAP Mask and Fit Pack

The **DreamEasy 2 Full Face CPAP Mask** provides a superior seal with its removable Comfort Cushion that conforms to the facial contours around the nose and mouth. Equipped with a 360° rotational elbow and standard 22 mm tubing connector allowing for freedom of movement when sleeping. The unique keyhole design forehead bar provides stability while eliminating pressure. Quick-release clips allow for easy removal of the mask without re-adjusting the headgear. Dishwasher safe for easy cleaning (top rack only).



FEATURES & BENEFITS

360° Rotational Elbow

Allows for more freedom of movement when sleeping with standard 22 mm tubing.

Eliminate Pressure

The unique keyhole design forehead bar provides stability while eliminating pressure.

Quick-Release

Clips allow for easy removal of the mask without re-adjusting the headgear.

Dishwasher Safe

for easy cleaning (top rack only).

Full Face Design

Perfect for patients whose therapy is compromised by mouth breathing.

Lightweight Design

Durable and lightweight designed face mask with hook-and-loop fastened headgear allows for maximum movement.

Comfort Cushion

Made with medical liquid silicone, the cushion is removable, odorless and provides an airtight seal.

CO₂ Washout

Elbow vent holes provide effective CO₂ washout.



Compass Health offers a 30-day, no questions asked, money-back guarantee on all Roscoe CPAP Masks, effective from the original date of purchase.

| Item # | Description | Size | HCPCS Code | Qty |
|------------|--|--------------------|------------|------|
| CPM-DEF2S | DreamEasy 2 Full Face CPAP Mask with Headgear | Small | A7030 | Each |
| CPM-DEF2M | DreamEasy 2 Full Face CPAP Mask with Headgear | Medium | A7030 | Each |
| CPM-DEF2L | DreamEasy 2 Full Face CPAP Mask with Headgear | Large | A7030 | Each |
| CPM-DEF2SK | DreamEasy 2 Full Face CPAP Mask with Headgear | All Sizes Fit Pack | A7030 | Each |
| X1DEF2HS01 | DreamEasy 2 Full Face Replacement Seal | Small | A7031 | Each |
| X1DEF2HS02 | DreamEasy 2 Full Face Replacement Seal | Medium | A7031 | Each |
| X1DEF2HS03 | DreamEasy 2 Full Face Replacement Seal | Large | A7031 | Each |
| X1DEF2HS04 | DreamEasy 2 Full Face Replacement Headgear, Gray | One Size | A7035 | Each |
| X1DEF2HS05 | DreamEasy 2 Full Face Replacement Headgear Clips | N/A | N/A | Pair |

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

Roscoe Medical DreamEasy™ Full Face Mask

The DreamEasy™ Full Face CPAP Mask provides a superior seal with its removable Comfort Cushion™ that conforms around the nose and facial contours for a comfortable night's rest. Equipped with a 360° rotational elbow and standard 22 mm tubing connector allowing for freedom of movement when sleeping. Comfortable forehead silicone pad diffuses any pressure. Quick-release clips allow for easy removal of the mask without re-adjusting the headgear. Dishwasher safe for easy cleaning. Replacement Comfort Cushion seals available.



PRODUCT FEATURES

- Full face design is perfect for patients whose therapy is compromised by mouth breathing
- Extremely lightweight and durable
- Headgear with Velcro® fasteners provide maximum versatility and movement
- 360° rotational elbow with standard 22 mm tubing connector
- The Comfort Cushion™ made with medical liquid silicone is completely odorless
- Comfort Cushion can be removed for easy cleaning
- 90 day mask and Comfort Cushion kits available
- Ergonomic design is airtight
- Vent hole provides effective CO₂ washout
- Latex-free
- Mask is dishwasher safe (top rack) for easy cleaning
- 90 day limited warranty

| ITEM # | DESCRIPTION | QTY. |
|-----------|--|------|
| CPM-DEFS | DreamEasy™ Full Face Mask with Headgear, Small | each |
| CPM-DEFM | DreamEasy™ Full Face Mask with Headgear, Medium | each |
| CPM-DEFL | DreamEasy™ Full Face Mask with Headgear, Large | each |
| 90108 | DreamEasy™ Full Face Mask Replacement Seal, Small | each |
| 90110 | DreamEasy™ Full Face Mask Replacement Seal, Medium | each |
| 90112 | DreamEasy™ Full Face Mask Replacement Seal, Large | each |
| 90395 | Replacement Headgear Clips | 2/pk |
| 90377 | Replacement Headgear | each |
| CPAP-SZGD | Nasal/Full Face Mask Sizing Guide | each |

BILLING TIPS (see page 19)

HCPCS CODES

A7030: Mask Interface,
Full Face Style

A7031: Replacement Full Face
Cushion

A7035: Headgear



Roscoe Medical® offers a 30-Day Money-Back Guarantee on all Roscoe® CPAP Masks, effective from the date of purchase.

If your patient is not 100% completely satisfied and it is within 30 days from your Roscoe purchase date, we will refund the purchase price to you, no questions asked.

Roscoe Medical® DreamEasy™ Nasal Mask

The DreamEasy™ Nasal CPAP Mask provides a superior seal with its removable Comfort Cushion™ that conforms around the nose and facial contours for a comfortable night's rest. Equipped with a 360° rotational elbow and standard 22 mm tubing connector allowing for freedom of movement when sleeping. Comfortable forehead silicone pad diffuses any pressure. Quick-release clips allow for easy removal of the mask without readjusting the headgear. Dishwasher safe for easy cleaning. Replacement Comfort Cushion seals available.



PRODUCT FEATURES

- Flexible and lightweight design allows for a comfortable, secure fit
- Headgear with Velcro® fasteners provide maximum versatility and movement
- 360° rotational elbow with standard 22 mm tubing connector
- The Comfort Cushion™ made with medical liquid silicone is completely odorless
- Comfort Cushion can be removed for easy cleaning
- Ergonomic design is airtight
- Vent hole provides effective CO₂ washout
- Latex free
- Mask is dishwasher safe (top rack) for easy cleaning
- 90-day limited warranty

| ITEM # | DESCRIPTION | QTY. |
|-----------|--|------|
| CPM-DENS | DreamEasy™ Nasal Mask with Headgear, Small | each |
| CPM-DENM | DreamEasy™ Nasal Mask with Headgear, Medium | each |
| CPM-DENL | DreamEasy™ Nasal Mask with Headgear, Large | each |
| CPM-DENSK | DreamEasy™ Nasal Mask with Headgear, Starter kit | each |
| CPAP-SZGD | Nasal/Full Face Mask Sizing Guide | each |
| 90091 | DreamEasy™ Replacement Nasal Seal, Small | each |
| 90093 | DreamEasy™ Replacement Nasal Seal, Medium | each |
| 90095 | DreamEasy™ Replacement Nasal Seal, Large | each |
| 90395 | Replacement Headgear Clips | 2/pk |
| 90377 | Replacement Headgear | each |

BILLING TIPS (see page 19)

HCPCS CODES

A7034: Mask Interface, Nasal or Cannula Style
A7032: Replacement Nasal Cushion
A7035: Headgear

Roscoe Medical® Zzz - Mask Full Face CPAP Mask

The Roscoe® Full Face CPAP Mask's one piece design is flexible, lightweight, and very easy to clean. The soft medical grade silicone mask allows for a comfortable fit and excellent seal without the need of a forehead support. Equipped with a 360 degree flexible rotational elbow and standard 22 mm tubing connector allowing for freedom of movement when sleeping.



PRODUCT FEATURES

- Made of high-quality, lightweight, clear silicone with a flexible frame
- Dual 360° swivel at the tubing attachment point
- Headgear with Velcro® fasteners provide maximum versatility
- Latex free
- Dishwasher safe (top rack) for easy cleaning
- Full face mask comes complete with headgear (black)
- 90-day limited warranty

| ITEM # | DESCRIPTION | QTY. |
|-----------|---|------|
| PB7800S | ZZZ-Mask Full Face Mask with Headgear, Small | each |
| PB7800M | ZZZ-Mask Full Face Mask with Headgear, Medium | each |
| PB7800L | ZZZ-Mask Full Face Mask with Headgear, Large | each |
| CPAP-SZGD | Nasal/Full Face Mask Sizing Guide | each |
| PB780S | ZZZ-Mask Full Face Seal, Small | each |
| PB780M | ZZZ-Mask Full Face Seal, Medium | each |
| PB780L | ZZZ-Mask Full Face Seal, Large | each |

BILLING TIPS (see page 19)

HCPCS CODES


A7030: Mask Interface, Full Face Style
A7035: Headgear
A7031: Replacement Full Face Cushion

Roscoe Medical® Shadow™ Nasal Pillows CPAP Mask

The Roscoe® Shadow™ Nasal Pillows CPAP Mask features a sleek design that conforms naturally to the contours of the face providing a comfortable fit with minimal facial contact. This permits for almost any sleep position without mask displacement. This innovative mask allows the user to anchor the tube at the top of the headgear or to leave it free to track. The anatomically shaped soft silicon pillows range in size from extra-small to large for a perfect fit. A complete CPAP solution right out of the box! **Replacement seals available.**

PRODUCT FEATURES

- Lightweight sleek design conforms to facial contours with minimal contact for maximum comfort and compliance
- Anatomically shaped soft silicone pillows provide maximum comfort with an excellent seal
- 360° rotational elbow design allows for tubing can be anchored at the top of the headgear or left free to track
- A built-in air vent disperses air quietly to prevent sleep disturbance
- Fit pack comes complete with 4 sets of pillows: extra-small, small, medium and large
- Replacement pillows available
- Latex free
- 90-day limited warranty



| ITEM # | DESCRIPTION | QTY. |
|--------|---|------|
| 70116 | Shadow™ Nasal Pillows Mask Fit Pack with 4 Sets of Pillows and Headgear | each |
| 90385 | Shadow™ Replacement Pillows, extra small | pair |
| 90387 | Shadow™ Replacement Pillows, small | pair |
| 90389 | Shadow™ Replacement Pillows, medium | pair |
| 90391 | Shadow™ Replacement Pillows, large | pair |

| BILLING TIPS (see page 19) |
|---|
| HCPCS CODES |
| A7034: Mask Interface, Nasal or Cannula Style |
| A7033: Replacement Pillows |
| A7035: Headgear |

Sleepnet® Ascend® CPAP Mask Series

The Sleepnet® Ascend® CPAP Mask Series features advanced cushion technology with the replaceable AIR°gel® cushion. For added convenience, patients can replace the cushion without having to replace the whole mask. The patented AIR°gel provides maximum comfort and minimum weight helping to reduce pressure points and skin irritation for better patient compliance. The flexible spacebar has a unique bendable design that provides the patient with superior stability without contacting the skin.



Ascend® Nasal Mask

Ascend® Full Face Mask

PRODUCT FEATURES

- Detachable cushion allows you to replace your cushion without replacing your entire mask
- Flexible spacebar is designed to provide superior stability without contacting the skin
- Active headgear connector moves when the patient moves to help maintain an effective seal
- Available in 3 sizes (small, medium and large) with simple and breathable quick-release EZ Fit headgear

- 360° rotating swivel for optimal tube positioning and easier mobility
- Patented AIR°gel® with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal
- Available in your choice of full face or nasal mask

| ITEM # | DESCRIPTION | QTY. |
|--------|---|------|
| 50825 | Sleepnet® Ascend® Full Face Mask System with Headgear, small, medium and large cushion included | each |
| 50207 | AIR°gel® Full Face Replacement Cushion, Small | each |
| 50208 | AIR°gel® Full Face Replacement Cushion, Medium | each |
| 50209 | AIR°gel® Full Face Replacement Cushion, Large | each |
| 50174 | Sleepnet® Ascend® Nasal Mask System with Headgear, small, medium and large cushion included | each |
| 50201 | AIR°gel® Nasal Mask Replacement Cushion, Small | each |
| 50202 | AIR°gel® Nasal Mask Replacement Cushion, Medium | each |
| 50203 | AIR°gel® Nasal Mask Replacement Cushion, Large | each |

BILLING TIPS (see page 19)

HCPCS CODES

- A7034: Mask Interface, Nasal or Cannula Style
- A7030: Mask Interface, Full Face Style
- A7031: Replacement Full Face Cushion
- A7032: Replacement Nasal Cushion
- A7035: Headgear

Sleepnet® MiniMe® 2 Pediatric Nasal Mask

Designed specifically for children, the MiniMe® 2 is available in two sizes for a greater range of comfort for children between ages two and twelve. Enhanced Custom Fit Technology™ enables the mask to be easily molded to a child's face maximizing comfort. The MiniMe 2 also features Advanced Cushion Technology™ (A.C.T.), Sleepnet's new ergonomic cushion technology delivers an effective seal yet minimizes pressure on small faces.



PRODUCT FEATURES

- Specifically engineered with children in mind, the Sleepnet® MiniMe® 2 provides unmatched comfort, a non-invasive and customized design, easy adjustment, and a pediatric size frame for the extra small faces we care so much about
- Touchless Spacebar™ eliminates any contact with the patient's forehead
- Custom Fit Technology™ means the shell of the mask can be gently molded to match your needs

- Designed specifically for children, the MiniMe® 2 is available in 2 sizes for a greater range of comfort for children ages 2–12
- Lightweight, breathable headgear
- Designed and built in the USA
- Does not contain any natural rubber latex

| ITEM # | DESCRIPTION | QTY. |
|--------|--|------|
| 60217 | MiniMe® 2 Pediatric Mask with Headgear, Small/Medium | each |
| 60218 | MiniMe® 2 Pediatric Mask with Headgear, Medium/Large | each |

BILLING TIPS (see page 19)

HCPCS CODES

- A7034: Mask Interface, Nasal or Cannula Style
- A7035: Headgear

Resvent iBreeze Auto CPAP & Bilevel Systems

Introducing the Resvent iBreeze Auto CPAP and Bilevel Systems with Heated Humidification, feature-rich PAP devices that were granted EUA (Emergency Use Authorization) by the FDA in 2020.* Some of the features include IPR (Intelligence Pressure Release), an algorithm during inspiration and expiration offering patients a more comfortable breathing experience; an auto humidifier which adjusts the heat based on room temperature and humidity to prevent rainout in the hose or mask; and finally, a 2-in-1 humidifier design which monitors the water level and when low, will signal a visible alert and automatically switch off the heat to prevent damage to the device. The iBreeze Systems come complete with SD Card, Carry Bag, 6' CPAP Tube, Tube Connector, Power Cord, and User Manual. *Wi-Fi equipped devices optional.*



SPECIFICATIONS

| | |
|-----------------------|---|
| H2O Therapy Pressure: | 4-20 cm (APAP) 4-30 cm (BPAP) |
| Ramp Time: | 0-60 minutes (5-minute increments) |
| Dimensions: | 9.25" w x 5.8" d x 4.9" h (device); 14.25" x 9.75" x 7" (shipping) |
| Weight: | 3.4 lbs (device); 6.2 lbs (shipping) |
| Water Capacity: | 290 mL |
| Humidifier Settings: | Auto / 0-8 (Default: 3) |
| Humidifier Output: | ≥10 mg H2O/L |
| Warranty: | Two Years Limited (device) |



FEATURES & BENEFITS

Large 3.5" LCD Display

Easily adjust settings and view last session to encourage user compliance.

Quiet Operation

with real time audible alert for tubing/mask leaks and automatic leakage/altitude compensation.

Intelligence Pressure Release

with 3 comfort settings and inspiration trigger for automatic startup.

Compliance Data

access compliance reporting via SD card or optional Wi-Fi.

Embedded Memory

stores previous night and 365 night user raw data.

iMatrix/ResAssist

Register your patients to maintain records and generate sleep reports on iMatrix desktop software or ResAssist cloud platform.

| ITEM | DESCRIPTION | HCPSCS | UOM | QTY |
|-------------------|--|--------|------|-----|
| IBREEZE20A | iBreeze Auto CPAP System with Heated Humidification | E0601 | Each | 1 |
| IBREEZE20A-W | iBreeze Auto CPAP System with WiFi and Heated Humidification | E0601 | Each | 1 |
| IBREEZE30STA | BiLevel Auto S/T System with WiFi and Heated Humidification | E0471 | Each | 1 |
| Replacement Parts | | | | |
| 05205-000001-00 | iBreeze Replacement 2-Piece Power Adapter with Brick | N/A | Each | 1 |
| 1001-000222-00 | iBreeze Replacement Humidifier Chamber | A7046 | Each | 1 |
| CPF-IBREEZE | iBreeze Replacement Disposable Filter | A7038 | Pack | 2 |
| 69999-000032-00 | iBreeze APAP Replacement Knob | N/A | Each | 1 |
| 69999-000033-00 | iBreeze BPAP Replacement Knob | N/A | Each | 1 |
| 10001-010043-00 | iBreeze Replacement Tube Adapter | N/A | Each | 1 |

Roscoe Medical® Private Label Program

With Roscoe Medical’s CPAP Private Label Program, you have the ability to build your brand awareness with CPAP replenishment supplies and improve your repeat patient business. Create a professional image and improve patient loyalty by private labeling your CPAP tubing, filters and customized tubing/filter kits. Include an insert with helpful replacement tips and related cash items, all branded to your company with your information. You will enjoy repeat supply orders and additional cash sales directly from your patients.



PROGRAM BENEFITS

- Low minimum order quantity of 50 each/pack with quick turnaround time
- Grow your repeat business for replenishment supplies
- Boost cash sales for CPAP-related accessories (Volt™ Battery Pack, pillows, wipes, REMZzzz's, hose wraps, etc...)
- Helps maintain profitability in the wake of Round 2 reimbursement cuts
- Provides the ability to brand replenishment supplies such as:
 - › CPAP Tubing
 - › Filters
 - › Tubing/Filter Kits
- Create a professional image for your business
- Include helpful replacement tips informing patients how often tubing and filters should be replaced



Sample Private Label




DreamEasy2 Full Face Mask & Easy-Flex6 Tube



Ascend Full Face Mask, Easy-Flex6 Tube, O2 Bleed-in Adapter & 25' Oxygen Tube

| SAMPLE PRIVATE LABEL KIT | | QTY |
|--------------------------|--|-----|
| CPM-DEF2SK | DreamEasy2 Full Face Mask | 1 |
| EASY-FLEX6 | Easy-Flex 6' CPAP Tubing | 1 |
| 50825 | Sleepnet Ascend Full Face Mask w/ EZFit Headgear | 1 |
| AG1642 | O2 Bleed-in Adapter | 1 |
| TUB-ROS25G | 25' Oxygen Tube Green | 1 |




iBreeze Filters For Resvent Units

| ITEM # | | UNIT | QTY | DIMENSIONS | OEM | AG | SUNSET | HCPCS CODE |
|-------------|---|---------|------|----------------------|-----|-----|--------|------------|
| CPF-IBREEZE |  | Resvent | 2/pk | 6.1" x 3.94" x 0.47" | N/A | N/A | N/A | A7038 |




Disposable Filters For Philips Respironics Units

| ITEM # | | UNIT | QTY | DIMENSIONS | OEM | AG | SUNSET | HCPCS CODE |
|------------|---|--|------|--------------------------------|---------|--------------|------------|------------|
| CPF-DSD2 |  | DreamStation | 2/pk | 2 1/4" x 1" | 1122518 | N/A | N/A | A7038 |
| CPF-PR1 |  | PR System One | 2/pk | 1 3/4" x 7/8" | 1063096 | AG1063096MED | CF1009NT-2 | A7038 |
| CPF-F30PK2 |  | Remstar Pro M Series, SleepEasy | 2/pk | 1 3/4" x 1 1/4" (includes tab) | 1029331 | AG1029331M | CF1007-2 | A7038 |
| CPF-F19 |  | Remstar Plus, Pro, BiPAP Pro 2 with Bi-Flex, Remstar Autosmart | 2/pk | 3 5/8" x 2" (includes tab) | 1005945 | AG5945 | CF1006 | A7038 |


Reusable Filters For Philips Respironics Units

| ITEM # | | UNIT | QTY | DIMENSIONS | OEM | AG | SUNSET | HCPCS CODE |
|------------|---|--|------|----------------------|---------|------------|-----------|------------|
| CPF-A11 | NO IMAGE AVAILABLE | Airsense 11 | 2/pk | 4.65" x 3.78" x 1.5" | 36851 | AG36851 | CF2107-2 | A7038 |
| CPF-DSR1 |  | Dream Station | 1/pk | 2 1/4" x 7/8" x 5/8" | 1122446 | A7039 | N/A | N/A |
| CPF-F31PK1 |  | PR System One, Remstar M Series, SleepEasy | 1/pk | 1 3/4" x 7/8" | 1029330 | AG1029330M | CF1007F-1 | A7039 |
| CPF-F28PK1 |  | Remstar Plus, Pro, BiPAP Pro 2 w/Bi-Flex, Remstar Autosmart, Legacy C-Flex | 1/pk | 3 5/8" x 1 5/8" | 1005964 | AG5964M | CF1006F | A7039 |



Filters For ResMed® Units

| ITEM # | | UNIT | QTY | DIMENSIONS | OEM | AG | SUNSET | HCPCS CODE |
|-------------|---|--|------|----------------------|-------|------------|----------|------------|
| CPF-A11 | NO IMAGE AVAILABLE | Airsense 11 | 2/pk | 4.65" x 3.78" x 1.5" | 36851 | AG36851 | CF2107-2 | A7038 |
| CPF-S9PK2 |  | S9 Series, AirSense 10, AirCurve 10 | 2/pk | 2 1/8" x 1 7/16" | 36851 | AG36851 | CF2107-2 | A7038 |
| CPF-S9HYPK2 |  | S9 Series Hypoallergenic, AirSense 10, AirCurve 10 | 2/pk | 2 1/8" x 1 7/16" | 36856 | AG36856 | CF2107HA | A7038 |
| CPF-S8PK2 |  | S8 Series | 2/pk | 1 5/16" x 1 3/8" | 33916 | AF33916MED | CF2105-2 | A7038 |


Filters for DeVilbiss® Units

| ITEM # | | UNIT | QTY | DIMENSIONS | OEM | AG | SUNSET | HCPCS CODE |
|-------------------|---|--|------|------------------|-----------|---------|--------|------------|
| ULTRA FINE FILTER | | | | | | | | |
| CPF-F32PK2 |  | IntelliPAP Standard, Standard Plus, Auto Adjust and Bi-level S | 2/pk | 1 5/8" x 1 3/16" | DV51D-603 | F603MED | CF5004 | A7038 |

Filters for 3B Medical Units

| ITEM # | | UNIT | QTY | DIMENSIONS | OEM | AG | SUNSET | HCPCS CODE |
|-----------|---|--------------|------|-----------------|------|-----|----------|------------|
| LG1020 |  | Luna/Luna II | 6/pk | 2" x 1 1/4" | Same | N/A | CF4003-2 | A7038 |
| CPF-3BPK2 |  | RESmart | 2/pk | 3 1/4" x 1 5/8" | N/A | N/A | CF8004 | A7038 |

Bacteria Filters

| ITEM # | | UNIT | QTY | DIMENSIONS | OEM | AG | SUNSET | HCPCS CODE |
|--------|---|----------|------|------------|------|------|--------|------------|
| AG7178 |  | Multiple | each | 22 mm | Same | Same | BF2200 | N/A |

A photograph of an older couple walking through a sun-dappled forest. The woman, on the left, has blonde hair and is wearing a blue t-shirt. The man, on the right, has grey hair and glasses, wearing a light blue polo shirt and khaki pants. They are both smiling and looking towards the right. The background consists of many trees with green leaves, and the ground is covered in fallen leaves and dappled sunlight.

CPAP TUBING & ACCESSORIES

Our Roscoe Medical® high quality CPAP Tubing is offered in a variety of styles. The tubing features a smooth interior which allows for excellent air flow and is compatible with virtually all CPAP, APAP and Bi-Level Systems. We also offer a wide range of headgear and accessories to meet all your patient's needs.

**Please note: All Roscoe Medical® CPAP tubing is intended for single-patient use.*

Easy-Flex™ Lightweight CPAP Tubing

The Roscoe Medical® Easy-Flex™ Lightweight CPAP Tubing has a soft, flexible feel providing patients with the high quality they expect. Equipped with 22 mm cuffs, our improved tubing is lighter and more flexible than comparable tubing. This tubing is dark gray in color and available in a variety of lengths.



#EASY-FLEX6

| ITEM # | DESCRIPTION | QTY. |
|--------------|--|------|
| EASY-FLEX6* | Easy-Flex™ Lightweight Tubing, 6 ft. Dark Gray with 22 mm Cuffs | each |
| EASY-FLEX8* | Easy-Flex™ Lightweight Tubing, 8 ft. Dark Gray with 22 mm Cuffs | each |
| EASY-FLEX10* | Easy-Flex™ Lightweight Tubing, 10 ft. Dark Gray with 22 mm Cuffs | each |
| EASY-FLEX18* | Easy-Flex™ Lightweight Tubing, 18 in. Dark Gray with 22 mm Cuffs | each |

BILLING TIPS (see page 19)

HCPCS CODES

A7037: CPAP Tubing

FlexiSlim™ CPAP Tubing

The Roscoe Medical® FlexiSlim™ CPAP Tubing is a small diameter air tube that is exceptionally light and flexible. FlexiSlim tubing works well with smaller lightweight masks as it is less likely to pull on the mask and break the seal.



PRODUCT FEATURES

- 6 ft. dark gray, small diameter air tube
- Ergonomic 22 mm cuffs for easy installation and removal from CPAP unit
- Exceptionally lightweight and flexible—40% lighter than standard tubing
- FlexiSlim™ tubing is 10% thinner than standard CPAP tubing
- Patient can freely move around in bed without creating mask leaks
- Tubing is designed to work with ResMed S9/S10 series and any devices with 15 mm tubing settings
- 90-day warranty

| ITEM # | DESCRIPTION | QTY. |
|------------|--|------|
| SLIM-FLEX6 | FlexiSlim™ CPAP Tubing, 6 ft. Dark Gray with 22 mm Cuffs | each |

BILLING TIPS (see page 19)

HCPCS CODES

A7037: CPAP Tubing

CPAP REPLACEMENT HEADGEAR & ACCESSORIES

Patients may need to replace their headgear occasionally as it may become stretched or lose their elasticity. For optimal performance, provide your patients your patients with a personalized replacement plan to ensure optimal performance and help keep them compliant.

CPAP Replacement Headgear



#ROS-FFHG



#ROS-NHG

BILLING TIPS (see page 19)

HCPCS CODES

A7035: Headgear

| ITEM # | DESCRIPTION | QTY. |
|------------|--|------|
| X1DEF2HS04 | DreamEasy 2 Full Face Replacement Headgear, Gray | each |
| X1DEF2HS05 | DreamEasy 2 Full Face Replacement Headgear Clips | pair |
| 90377 | Roscoe Medical® Headgear, for use with Roscoe® DreamEasy™ Full Face and Nasal CPAP Masks, Black | each |
| ROS-NHG | Universal Nasal Mask Headgear, Black | each |
| ROS-FFHG | Universal Full Face Mask Headgear, Black | each |
| 90395 | Roscoe Medical® Replacement Headgear Clips, for use with Roscoe® DreamEasy™ Full Face and Nasal CPAP Masks | 2/pk |

CPAP Replacement Chin Straps



#ROS-T09M

| ITEM # | DESCRIPTION | QTY. |
|--------------|--|------|
| ROS-T09S | 3-Point Chin Strap, Small (compare to Tiara/Cardinal #TMS-09S), Ruby Red | each |
| ROS-T09M | 3-Point Chin Strap, Medium (compare to Tiara/Cardinal #TMS-09M), Ruby Red | each |
| ROS-T09L | 3-Point Chin Strap, Large (compare to Tiara/Cardinal #TMS-09L), Ruby Red | each |
| ROS-T09ADJ | 3-Point Adjustable Chin Strap,with adjustable 2 in. extension (compare to Tiara/Cardinal #TMS-09ADJ), Ruby Red | each |
| ROS-T09ADJXL | 3-Point Adjustable Chin Strap, Extra Large, Ruby Red | each |
| AG1012911 | Premium Chin Strap (compare to Respirationics #1012911), Ruby Red | each |
| AC302175B | Philips Respirationics-Style Chin Strap, White | each |
| AG302425 | Philips Respirationics-Style Deluxe Chin Strap (compare to #302425) | each |
| AC133318 | Puritan Bennett™-Style Chin Strap (compare to Puritan Bennett™ #133318), Ruby Red | each |

CPAP Mask Wipes

These 100% pure cotton wipes are specially formulated to naturally clean and protect your mask without the use of harmful chemicals. CPAP Mask Wipes effectively remove dirt, grease, oils and other residue from the patient's mask keeping it fresh and clean. Wipes are made from 100% natural ingredients and contain no strong or unpleasant fragrances.

PRODUCT FEATURES

- Contain only natural ingredients derived from coconut and other plant extracts
- Suitable for daily cleaning of CPAP masks, tubing and accessories
- Gentle on hands and face
- 100% pure cotton
- Alcohol and latex free
- Biodegradable



FGC15600

FGC13700

| ITEM # | DESCRIPTION | QTY. |
|---------------|--|------|
| FGC13700 0000 | Carex CPAP Mask Wipes, Unscented (62 wipes per tub) | each |
| FGC13800 0000 | Carex CPAP Mask Wipes, Citrus Scent (62 wipes per tub) | each |
| FGC15600 0000 | Carex CPAP Mask Wipes Travel Pack, Unscented (10 individually wrapped wipes per box) | each |

CPAP Tubing Hose Wraps

Roscoe Medical® Hose Wraps eliminate the most annoying parts of therapy—the cold medical tubing that connects the user's system to the mask or interface. By adding our Roscoe® Hose Wraps to your patient's system, you can help eliminate condensation and rainout that often disrupts sleep and therapy.

PRODUCT FEATURES

- Lightweight and easy to travel with
- Ultra soft fleece material
- Easy to install and very patient friendly
- Fits any 6 ft. 22 mm CPAP tube
- Comes standard with easy grip closure
- Assists in controlling rainout
- Machine washable for easy cleaning
- Available in 4 colors: burgundy, navy, olive and pink



#HW-NAVY

| ITEM # | DESCRIPTION | QTY. |
|---------|-------------------------|------|
| HW-NAVY | Roscoe® Hose Wrap, Navy | each |

CPAP Adaptors, Port and Humidification Chambers



AG1642



HC325S



HC385S



1001-00222-00

| ITEM # | DESCRIPTION | QTY. |
|---------------|---|------|
| 1001-00222-00 | iBreeze Replacement Humidifier | each |
| AG1642 | Oxygen Bleed-In Port | each |
| HC325S | Humidification Chamber, for HC150 CPAP Systems, Single-Patient Reusable | each |
| HC385S | Standard Humidification Chamber for SleepStyle™ 200 CPAP Series | each |

Carex CPAP Pillow

The Carex CPAP Pillow can help reduce CPAP mask leaks, reduces mask pressure and helps increase comfort for a good night's rest. The larger cutouts on either side allow the user to sleep on either side without worrying about leaks.

PRODUCT FEATURES

- Use with Full Face, Nasal or Direct Nasal interfaces
- Freedom to move without displacing mask
- Provides neck and head support
- Removable cover for easy cleaning




| ITEM # | DESCRIPTION | QTY. |
|---------------|--|------|
| FGP11300 0000 | Carex CPAP Pillow with Cover (box packaging) | each |

Contour CPAP Sleep Pillow

This specially designed pillow is for CPAP mask users and is used to improve comfort for all styles and brands of CPAP masks.

PRODUCT FEATURES

- Patented hollow areas alleviate mask pressure, improves comfort and helps prevent mask leaks
- Dual cutout design accommodates the mask, allowing for comfortable left- or right-side sleeping without compromising CPAP compliance
- Contoured shape supports head and neck improving spine alignment
- Improves airway alignment
- High-quality medical grade foam is cushiony yet holds its shape to prevent edge collapse, while maintaining the hollowed-out mask area
- Overall design helps CPAP users attain a comfortable, restorative sleep night after night



| ITEM # | DESCRIPTION | QTY. |
|--------|---------------------------------------|------|
| 70014 | Contour CPAP Pillow with Velour Cover | each |

Contour CPAPMax Pillow

The CPAPMax Pillow can help reduce CPAP mask leaks, reduces mask pressure and helps increase comfort for a good night's rest. The larger cutouts on either side allow the user to sleep on either side without worrying about leaks.

PRODUCT FEATURES

- Reduces mask interference, shifting and leaks
- Adjustable pillow height switched from regular to extra thick
- Minimizes facial pressure, sores and bruising
- Improve sleep comfort and CPAP compliance
- Ventilated foam is cool, breathes and reduces perspiration
- Ergonomic center cradles head and improves airway alignment
- Front edge matches the natural curve of neck and shoulders
- Pressure-free zones prevent interference between mask and pillow
- CoolTech 3D mesh cover
- 2-in-1 reversible memory foam or fiber top



| ITEM # | DESCRIPTION | QTY. |
|----------|------------------------|------|
| PLW-CPMX | Contour CPAPMax Pillow | each |

CPAP HCPCS Codes

| CODE | PRODUCT | REPLACEMENT SCHEDULE |
|-------|---|----------------------|
| A4604 | Tubing with integrated heating element for use with positive airway pressure (PAP) device | 1 per 3 months |
| A7027 | Combination oral/nasal mask, used with continuous positive airway pressure (CPAP) device, each | 1 per 3 months |
| A7028 | Oral cushion for combination oral/nasal mask, replacement only, each | 2 per month |
| A7029 | Nasal pillows for combination oral/nasal mask, replacement only, pair | 2 per month |
| A7030 | Full face mask used with PAP device, each | 1 per 3 months |
| A7031 | Face mask interface, replacement for full face mask, each | 1 per month |
| A7032 | Replacement cushion for nasal application device, each | 2 per month |
| A7033 | Replacement pillows for nasal application device, pair | 2 pairs per month |
| A7034 | Nasal mask or cannula type application device, used with PAP device | 1 per 3 months |
| A7035 | Headgear used with PAP device | 1 per 6 months |
| A7036 | Chin strap used with PAP device | 1 per 6 months |
| A7037 | Tubing used with PAP device | 1 per 3 months |
| A7038 | Filter, disposable, used with PAP device | 2 per month |
| A7039 | Filter, non-disposable, used with PAP device | 1 per 6 months |
| A7044 | Oral interface used with PAP device | Contact payor |
| A7045 | Exhalation port with or without swivel used with accessories for PAP devices, replacement only | Contact payor |
| A7046 | Replacement water chamber for humidification, used with PAP device | 1 per 6 months |
| A9279 | Monitoring feature/device stand-alone or integrated, any type, includes all accessories, components and electronics, not otherwise classified | Contact payor |
| E0470 | Respiratory assist device (RAD), bi-level pressure capability, without backup rate feature, used with noninvasive interface | Monthly rental rate |
| E0471 | Respiratory assist device, bi-level pressure capability, with backup rate feature, used with noninvasive interface | Monthly rental rate |
| E0561 | Humidifier, non-heated, used with PAP device | Contact payor |
| E0562 | CPAP heated humidifier | Contact payor |
| E0601 | CPAP device | Monthly rental rate |

*HCPCS Disclaimer:

HCPCS codes provided by Compass Health® are intended as general guidelines. Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended nor should it be considered billing advice. Compass Health® does not guarantee coverage or reimbursement on any product.



6753 Engle Road • Middleburg Hts., OH 44130 Toll Free: (800) 376-7263 • Fax: 440-572-4261

www.compasshealthbrands.com

rev_02_2023