

# OXYGEN CATALOG



# OXYGEN CONCENTRATORS

truAIRE-5  
by Roscoe Medical



## Stationary 5 Liter Oxygen Concentrator

Introducing the truAIRE-5 by Roscoe Medical, Compass Health Brands' first 5-liter stationary concentrator. The truAIRE-5 is a compact, low maintenance solution with easy, out-of-box set up and ultra-quiet operation. Designed to operate 24/7, the truAIRE-5 is the dependable 5-liter oxygen concentrator choice for your oxygen patients.

### Product Features:

- Compact and easy to transport
- Ultra-quiet operation
- Durable, double-wobble compressor designed for 24/7 use
- Low maintenance
- Easy removal and replacement of cabinet and compressor intake filters
- 5-Year Limited Warranty
- Includes Humidifier Adaptor

### Technical Specifications:

Dimensions:	13"L x 10.25"W x 21.25"H
Weight:	35.5 lbs
Liter Flow:	0.5 -5 LPM
Maximum Outlet Pressure:	5.5 psi $\pm$ 0.7 psi (35 kPa $\pm$ 5 kPa)
Oxygen Concentration:	93% ( $\pm$ 3%) at 0.5~5 LPM after warmup period of 30 minutes
Mode of Operation:	Continuous Flow
Safety System:	Low Pressure Alarm Low Oxygen Concentration Alarm High Pressure/High Temperature Shutdown Current Overload/Line Surge Shutdown
Electrical Requirements:	AC120V / 60 Hz
Rated Current Input:	3.5A
Sound Level:	$\leq$ 50.5 dBA
Warranty:	5 Year Limited Warranty

Item #	Description	Qty
O2C5L	truAIRE-5 Stationary 5 Liter Oxygen Concentrator	each



# Zen-O lite™ with RRT

## Portable Oxygen Concentrator

The Zen-O lite™ Portable Oxygen Concentrator offers reliability, flexibility and supports active lifestyles of patients requiring oxygen therapy. Zen-O lite™ features breath sensitivity and Rate Responsive Therapy™ (RRT) that automatically adjusts to the patient's demand.

### RATE RESPONSIVE THERAPY™ (RRT)

Zen-O lite™ delivers oxygen on demand which is triggered by the patient's inhalation. Whether a patient is climbing the stairs or taking a walk, Zen-O lite™ will automatically respond to the breath rate by adjusting the amount of oxygen needed to ensure the patient remains saturated.

### Product Features:

- Pulse Mode, up to 1050 ML
- Simple and Easy to Use
- Eco Mode
- Auto Mode
- Lightweight
- Easily replaceable sieve beds
- Breath detection indicator
- Approved for air travel
- Durable and reliable

### Technical Specifications:

Dimensions:	9.8" L x 3.8" W x 9.25" H
Weight:	5.5 lbs
Power Requirements:	AC Adaptor: 100-240V AC (+/- 10%) 50-60 Hz in, 24V DC, 5.0A out
	DC Adaptor: 11.5 - 16V DC in, 24V, 5.0A out
Purity:	87% -96% at all settings
Maximum oxygen discharge pressure:	20.5 psi
Inspiratory trigger sensitivity:	-0.12cm/H <sub>2</sub> O
Humidity Range:	5% to 93% ± 2% non-condensing
Operating Altitude	0m to 4,000m (0' to 13,000')
Operation Temperature:	41°F (5°C) and 104°F (40°C)
Storage Temperature:	-4°F (-20°C) and 140°F (60°C)
Setting:	Adjustable in 0.5 increments from 1.0 to 5.0
Noise Level:	37 dB(A)*
Battery Duration	Up to 4 hours*

\*At setting 2

Setting	Approx. Battery Life 1 Battery (Hours:Minutes)**
Pulse 1.0	4:00
Pulse 2.0	4:00
Pulse 3.0	3:00
Pulse 4.0	2:15
Pulse 5.0	2:00

\*\*Battery run times above are estimates based on 18 bpm

Item #	Description	Item #	Description
RS-00608-X-S	Zen-O lite™ portable concentrator with one 8 cell battery	RS-00608-X-D	Zen-O lite™ portable concentrator with two 8 cell batteries
RS-00601	Zen-O lite™ rechargeable battery	RS-00604	Zen-O lite™ AC power supply with US cord
RS-00605	Zen-O lite™ DC power supply	RS-00606	Zen-O lite™ carry bag
RS-00616	Zen-O lite™ replacement sieve beds	RS-00617	Zen-O lite™ cannula wrench
RS-00619	Zen-O lite™ rucksack	RS-00512	Zen-O lite™ cannula filter pk of 10
RS-00515	Zen-O lite™ single external battery charger US	RS-00523	Accessories bag
RS-00530	Zen-O lite™ dual external battery charger US		



## Medical Gas

Roscoe Medical is helping you give your oxygen patients the active lifestyle they want with our complete inventory of cylinders, conserving devices, regulators and cylinder supplies.



RMI-04ST



RMI-08MINI



RMI-15H

## Regulators and Conserving Devices

### Roscoe Medical Regulator

The Roscoe Medical Click-Style Regulator is constructed with an aluminum body and all brass core, providing patients with the ultimate in portability and safety. The Roscoe Regulator offers HME dealers a reliable, high-quality regulator at an affordable price. Liter flow rates for 0-8 regulators include settings of: 0, .12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6 and 8. Liter flow rates for 0-15 regulators include settings of: 0, .12, .5, 1, 1.5, 2, 3, 4, 6, 8, 10 and 15. Two year warranty. HCPCS Code: E1353\*

Item #	Flow Rate	Inlet Conn.	Flow Rates	Outlet	Qty
RMI-04ST	1/32 to 4	870	Pediatric, Hose Barb (Liter flow rates include: 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1 1/2, 2, 3, 4)	Hose Barb/ DISS 1240	each
RMI-04H	1/32 to 4	540	Pediatric, DISS 1240 (Liter flow rates include: 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1 1/2, 2, 3, 4)	Hose Barb/ DISS 1240	each
RMI-08MINI	0 to 8	870	Mini, Hose Barb Liter Flow Rates: 0, .12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, and 8	Hose Barb/ DISS 1240	each
RMI-15MINI	0 to 15	870	Mini, Hose Barb Liter Flow Rates: 0, .12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8, 10 and 15	Hose Barb/ DISS 1240	each
RMI-15H	0 to 15	540	Standard, DISS 1240 Liter Flow Rates: 0, .12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8, 10 and 15	Hose Barb/ DISS 1240	each



198705RM

## Roscoe EasyPulse 5 Conserving Device

The Roscoe® EasyPulse 5 is a small, lightweight oxygen conserving device, weighing in at just over half a pound. The single dial control and ergonomically designed T-Handle/ Hand Tight Nut make setting and connecting the EasyPulse 5 easy for any patient or caregiver. The EasyPulse 5 has 5 pulse settings and a 2 LPM continuous flow setting. With a conserving ratio of up to 5.7:1, the EasyPulse 5 reduces deliveries of cylinders to patients and its proven pneumatic reliability, robust design and lack of batteries helps to reduce overall operating costs. Two year warranty.

The EasyPulse 5 utilizes *Controlled Minute Volume (CMV)* Oxygen Conserving Technology. The patented CMV Technology delivers a pulse algorithm wave form, superior trigger sensitivity and fast delivery time. All this is accomplished pneumatically with a single lumen cannula.

Item #	Liter Flow	Inlet Conn.	Description	Qty
198705RM	1 to 5	CGA-870	Roscoe® EasyPulse 5 Oxygen Conserving Device	each



198712

## EasyPulse 5 & 5+6 Conserving Devices

The Precision Medical EasyPulse 5+6 has 6 continuous flow settings, to complement the 1-5 pulse settings. Even with the added continuous flow settings, the EasyPulse 5+6 weighs in at a mere 9.9 oz and is only 4.75 inches long, making it still the smallest and lightest conserving regulator available. In the Pulse Mode, the EasyPulse 5+6 delivers a front loaded tapered wave form. All the oxygen is delivered in the first third of the inspiratory phase. With the added security of more continuous flow options, you have the most versatile oxygen conserving device available.

Item #	Liter Flow	Inlet Conn.	Description	Qty
198712	1 to 6	870	Precision Medical EasyPulse 5+6 Oxygen Conserving Device	each
198705	1 to 5	870	Precision Medical EasyPulse 5 Oxygen Conserving Device	each



## Private Labeling?


Private labeling options are available for cylinder orders with a choice of stamping or silk screen. Pallet minimums are required. For further details, call your Compass Sales Representative at 1-800-3-ROSCOE (376-7263).

## Cylinders

Compass Health offers empty aluminum oxygen cylinders in standard sizes. You can order units in packs of six or in competitively priced pallet quantities.

- All cylinders are available in both CGA-870 wrench and toggle valves (excluding the CGA-540 configured sizes)
- Cylinder finish is easy to clean and maintain
- Private labeling options available

Note: Cylinder size shown is not proportionally correct



Cylinder Size	Item #	Dimensions	Oxygen Capacity (Liters)	Valve	Valve Type	Qty
ML6	PX-8705-1T	4 1/2" x 8"	170 Liters	870	Toggle	6/pack
ML6	PX-8705-1TB	4 1/2" x 8"	170 Liters	870	Toggle	81/pallet
ML6	PX-8705-1W	4 1/2" x 8"	170 Liters	870	Wrench	6/pack
ML6	PX-8705-1WB	4 1/2" x 8"	170 Liters	870	Wrench	81/pallet
M6	PX-8700-1T	11 3/4" x 3"	164 Liters	870	Toggle	6/pack
M6	PX-8700-1TB	11 3/4" x 3"	164 Liters	870	Toggle	144/pallet
M6	PX-8700-1W	11 3/4" x 3"	164 Liters	870	Wrench	6/pallet
M6	PX-8700-1WB	11 3/4" x 3"	164 Liters	870	Wrench	144/pallet
M9	PX-8701-1T	10 1/2" x 4 1/2"	246 Liters	870	Toggle	6/pack
M9	PX-8701-1TB	10 1/2" x 4 1/2"	246 Liters	870	Toggle	81/pallet
M9	PX-8701-1W	10 1/2" x 4 1/2"	246 Liters	870	Wrench	6/pack
M9	PX-8701-1WB	10 1/2" x 4 1/2"	246 Liters	870	Wrench	81/pallet
D	PX-8702-1T	16" x 4 1/2"	425 Liters	870	Toggle	6/pack
D	PX-8702-1TB	16" x 4 1/2"	425 Liters	870	Toggle	81/pallet
D	PX-8702-1W	16" x 4 1/2"	425 Liters	870	Wrench	6/pack
D	PX-8702-1WB	16" x 4 1/2"	425 Liters	870	Wrench	81/pallet
E	PX-8703-1T	25" x 4 1/2"	679 Liters	870	Toggle	6/pack
E	PX-8703-1TB	25" x 4 1/2"	679 Liters	870	Toggle	81/pallet
E	PX-8703-1W	25" x 4 1/2"	679 Liters	870	Wrench	6/pack
E	PX-8703-1WB	25" x 4 1/2"	679 Liters	870	Wrench	81/pallet
M60	PX-5405-1B	23" x 7.3"	1738 Liters	540	Spigot	25/pallet
M90	PX-5407-1B	32.7" x 7.25"	2,549 Liters	540	Spigot	25/pallet
M	PX-5406-1B	36" x 8"	3455 Liters	540	Spigot	25/pallet





## Medical Gas Home Filling Cylinders

Compass Health offers empty aluminum oxygen cylinders with Home Filling Units in four sizes packaged to meet your requirements. Available in single pack or 6 pack quantities.

- CGA-870 toggle valve available
- Brushed aluminum with clear coat finish
- Ships empty
- For use with Home Filling units
- Free freight on 6 packs for orders shipped within the Continental United States only

Item #	Description	UOM	Qty
HFML6-EA	HomeFill-Style ML6 Cylinder with Toggle Valve and Transfill Port	Each	1
HFC-EA	HomeFill-Style M9/C Cylinder with Toggle Valve and Transfill Port	Each	1
HFD-EA	HomeFill-Style D Cylinder with Toggle Valve and Transfill Port	Each	1
HFE-EA	HomeFill-Style E Cylinder with Toggle Valve and Transfill Port	Each	1

## Roscoe Adjustable Oxygen Carts

Roscoe Adjustable Oxygen Carts are constructed of durable, high-strength steel and can accommodate either "D" or "E" size cylinders. Handle is height adjustable to accommodate most individuals. Large, easy-roll wheels allow for safe maneuvering. Includes gurney/mask hook, nonslip hand grip and one piece thumb screw. Available in silver vein finish or with chrome plating. Limited lifetime warranty. HCPCS Code: E1355\*



E-CARTPC



E2-CARTPC



E-FOOTPK



E-SCREW10


Item #	Description	Height Adjustable	Cylinder Capacity	Weight	Qty
E-CARTPC	D or E Cylinder Cart, Silver Vein Finish	38 ½" - 42 ½"	1	5 lb.	4/cs
E2-CARTPC	D or E Cylinder Cart, Silver Vein Finish	38 ½" - 42 ½"	2	9 lb.	2/cs
E-FOOTPK	Rubber Feet for E-CART, E-CARTPC, E2-CARTPC				2/pk
E-SCREW10	Thumb Screw for E-CART				10/pk

## Medium-Duty Cylinder Racks and Stands

Compass Health's Cylinder Racks and Stands feature an all-welded frame construction to ensure durability and strength. Individual spacing, positive cylinder placement and weight distribution promote safe storage. Equipped with bolt down bars and feet, our stands are available in a green powder-coated finish.



CR-M6PC6



CR-EPC12




CR-MST2

Cylinder Size	Item #	Description	Cylinder Capacity	Rows	Dimensions (H x D x W, in.)	Weight	Qty
M6	CR-M6PC6	M6 Cylinder Rack, with Bolt Down Bars, Green Powder-Coat	6	2 x 3	13 x 8 ½ x 9 ½	5 lb.	each
M6	CR-M6PC12	M6 Cylinder Rack, with Bolt Down Bars, Green Powder-Coat	12	3 x 4	17 ¾ x 13 x 9 ¾	9 lb.	each
C/D/E	CR-EPC6	C/D/E Cylinder Rack, with Bolt Down Bars, Green Powder-Coat	6	2 x 3	16 ½ x 11 x 9 ¾	8 lb.	each
C/D/E	CR-EPC12	C/D/E Cylinder Rack, Bolt Down Bars, Green Powder-Coat	12	3 x 4	21 ¾ x 17 ¾ x 10 ¼	14 lb.	each
M60/ M90,MM	CR-MST2	M60, M90 and MM Cylinder Stand, Green Powder-Coat	1	n/a	14 ½ x 13 ¾ x 11 ¼	5 lb.	2/cs
D/E	150-0110	D/E Wheelchair Cylinder Carrier	1	n/a			4/cs
C/D/E	150-0265	Composite Cylinder Rack	8	4 x 2			each
M6	150-0275	Composite Cylinder Rack	12	4 x 3			each

## Cylinder Rack Parts and Accessories

Item #	Description	Qty
90014	Replacement Bolt for Roscoe CR-MST Cylinder Stand	each
150-0025	Vehicle Cylinder Restraint - Medium	each
X1JQMST01	Replacement Hand Nut for CR-MST	each



150-0025



## Heavy-Duty Cylinder Racks and Stands

These cylinder stands have been designed for use where there is a need for multiple cylinder storage and mobility is not a factor. The all-welded steel construction, combined with the tubular frame ensures greater durability, strength and longer life of the unit. Models 6060-LTL, 6120-LTL and 6240-LTL have a hinged lockable top to help prevent theft, misuse or mixing of D or E cylinders while in storage. Powder-coated finish.



150-0330



150-0320



150-0270



150-0290



150-0300



150-0250



610F-E

Cylinder Size	Item #	Description	Cylinder Capacity	Rows	Dimensions (H x D x W, in.)	Weight	Qty
M6	150-0330	M6 Cylinder Rack, with Bolt Down Bars, Green	15	3 x 4	19 ½ x 12 x 18	16 lb.	each
M6	150-0320	M6 Cylinder Rack, with Bolt Down Bars, Green	30	5 x 6	16 x 23 x 19	30 lb.	each
C/D/E	150-0270	C/D/E Cylinder Rack, with Bolt Down Bars, Green	12	3 x 4	19 ½ x 15 x 22	19 lb.	each
C/D/E	150-0290	C/D/E Cylinder Rack, with Bolt Down Bars, Green	24	4 x 6	19 ½ x 20 x 32	34 lb.	each
C/D/E	150-0300	C/D/E Cylinder Rack, with Bolt Down Bars, Green	40	5 x 8	20 x 40 x 25	40 lb.	each
M60/M/H	150-0250	M60/M/H Cylinder Stand, Black	1	n/a	11 x 18 ¼ x 19	16 lb.	each
M/M60/H	610F-E	M/M60/H Cylinder Stand, Holds 7" To 9 ½" Cylinder, Green	1	n/a	12 x 18 x 18	14 lb.	each
D/E	150-153	Horizontal Layered Rack	25	n/a	11 x 18 ¼ x 19	16 lb.	each
D/E	150-0153-DOT	Horizontal Layered Rack with Door	25	n/a	11 x 18 ¼ x 19	16 lb.	each

## Multiple Wheeled Cylinder Carts

Warehouse and Delivery Carts designed to withstand everyday use in high volume transport facilities. Equipped to transport a variety of cylinder sizes. Powder-coated finish.



150-0350



150-0140



150-0175



150-0235



150-0150



150-0158

Cylinder Size	Item #	Description	Cylinder Capacity	Rows	Weight	Qty
M4/M6/B	150-0350	Four 4" x 1 1/2" Casters (2 locking swivel, 2 rigid)	30	5 x 6	35 lb.	each
C/D/E	150-0140	Warehouse Cart with Four Caster Wheels	6	2 x 3		each
C/D/E	150-0170	Warehouse Cart with Four Caster Wheels	12	4 x 3		each
C/D/E	150-0175	Warehouse Cart with Two Caster Wheels and Two 10" Wheels	12	4 x 3		each
C/D/E	150-210	Warehouse Cart with Four Caster Wheels	24	6 x 4		each
C/D/E	150-0220	Warehouse Cart with Four Caster Wheels	28	7 x 4		each
MM/H	150-0235	Single Large Cylinder Transport Cart with Stair Slide	1	N/A		each
D/E	150-0150	Multi-Cylinder Delivery Cart, Horizontally Layered	68			each
D/E	150-0150-DOT	Multi-Cylinder Delivery Cart with Door & Top	68			each
D/E	150-0152	Horizontal Layered Cart	25			each
D/E	150-0152-DOT	Horizontal Layered Cart with Door	25			each
M6		Layered Cylinder Cart	84			each
C/D/E	150-0158	Vertical Transport Cart	40			each

## Medical Gas Accessories



1109-0014



PV-TG



PV-WR



ET-CW



2544-2



47012



47012



40150



SMK-SGN



1109-0005



150-2005PK10

Item #	Description	Qty
Cylinder Accessories		
PV-TG	Post valve, toggle-style, CGA-870	each
PV-WR	Post valve, wrench style, CGA-870	each
ET-CW	Cylinder wrench for wrench valve	25/pk
B-F66082	Cylinder wrench for H tank	10/bx
2544-2	Regulator yoke washer with rubber ring	25/pk
47012	Medical post valve locking band, white	each
40150	“EMPTY / IN USE / FULL” medical post valve tag	each
SMK-SGN	“No Smoking, Oxygen in Use” sign, non-adhesive, printed in red	100/pk
1109-0005	Cylinder labels (OXYGEN USP, UN1072), blank individual sheets,	each
Gauges		
1109-0014	Wet Band Post Valve Seals	1000/bucket
MPV-4046	Replacement Toggle and Pin	each
180-5005	Liter Meter, 1-8 LPM Pen-Style	each
180-5006	Liter Meter, 2-15 LPM Pen-Style	each
150-2005PK10	Disposable M6 Cylinder Tote, Holds 6	10/pk

## CYLINDER BAGS



37



31



24N



48N

Brand	Item #	Description	Qty
Zephyr Deluxe Black			
Zephyr	35	M6 Shoulder Bag	each
Zephyr	39	M7 Cylinder Shoulder Bag	each
Zephyr	33	C Cylinder Shoulder Bag	each
Zephyr	31	D Cylinder Shoulder Bag	each
Zephyr	29	M6/C M9 Cylinder Backpack	each
Airlift Premium Navy			
Airlift	44N	D/E Carrier For Wheelchair	each
Airlift	46N	Fanny Pack/Shoulder Bag For M6, M9/C	each
Airlift	48N	Fanny Pack/Shoulder Bag For ML6, M7	each
Airlift	49N	Camera Style Horizontal Shoulder	each
Airlift	24N	M6, M9/C Backpack/Shoulder	each
Airlift	26N	D Cylinder Backpack/Shoulder	each





## Oxygen Supply

Roscoe Medical's oxygen delivery products are designed with the highest quality materials for a comfortable patient fit, but at a lower price range than comparable product lines. Roscoe's Soft Nasal Cannulas feature an over-the-ear design that provide a high level of stability and patient comfort. From cannulas, tubing, humidifiers and more, Roscoe Oxygen Supply Products will satisfy all your patient needs.

## Nasal Cannulas

### Roscoe Medical Nasal Cannulas

Roscoe Medical Nasal Cannulas are comfortable and economical. Made of extra-soft disposable materials that are not made with natural latex rubber, Roscoe Cannulas are crush resistant and fit comfortably over the ears.



Item #	Description	Qty
CAN-PEDROS	Roscoe Pediatric Curved Soft Nasal Cannula with 7' Crush-Resistant Tubing	50/cs
CAN-ROS	Roscoe Adult Curved Soft Nasal Cannula without Supply Tubing	50/cs
CAN-ROS4	Roscoe Adult Curved Soft Nasal Cannula with 4' Crush-Resistant Tubing	50/cs
CAN-ROS7	Roscoe Adult Curved Soft Nasal Cannula with 7' Crush-Resistant Tubing	50/cs
CAN-ROS25	Roscoe Adult Curved Soft Nasal Cannula with 25' Crush-Resistant Tubing	25/cs

### Salter Nasal Cannulas

Salter's nasal cannulas feature a unique design that is anatomically curved to fit the upper lip with curved, tapered nasal prongs. It eliminates the flap or ridge which can cause irritation. Over-the-ear style ensures secure positioning, allowing complete freedom of movement. Demand cannula design provides separate independent tubes for oxygen delivery and sensing inspiratory effort.



1600



1615

Item #	Description	Qty
1600	Adult Cannula with 7' Three-Channel Safety Supply Tubing	50/cs
1600-4	Adult Cannula with 4' Three-Channel Safety Supply Tubing	50/cs
1600-14	Salter Cannula with 14' 3-Channel Supply Tubing	50/cs
1600-25	Adult Cannula with 25' Three-Channel Safety Supply Tubing	25/cs
1602	Pediatric Cannula with 7' Three-Channel Safety Supply Tubing (For Child One Year Of Age Or Older)	50/cs
1606	Adult Cannula without Supply Tubing	50/cs
1615	Intermediate Infant Cannula with 7' Three-Channel Safety Supply Tubing (For Child 3-12 Months Of Age)	50/cs
1616	Adult "Micro" (Low-Flow) Cannula with 7' Three-Channel Safety Supply Tubing	50/cs
1600HF-7-25	Adult High Flow Cannula with 7' Three-Channel Safety Supply Tubing	25/cs
16SOFT-4-50	Salter Soft Cannula with 4' 3-Channel Supply Tubing	50/cs
16SOFT-7-50	Salter Soft Cannula with 7' 3-Channel Supply Tubing	50/cs

## Salter Nasal Cannulas (Cont'd)

Item #	Description	Qty
1600HF-14-10	Salter High Flow Cannula with 14' 3-Channel Supply Tubing	10/cs
1600HF-25	Salter High Flow Cannula with 25' 3-Channel Supply Tubing	10/cs
1600HF-4-25	Salter High Flow Cannula with 4' 3-Channel Supply Tubing	25/cs
1600HF-50	Salter High Flow Cannula with 50' 3-Channel Supply Tubing	10/cs
1600HF-7-25	Salter High Flow Cannula with 7' 3-Channel Supply Tubing	25/cs
1056	Salter Cannula with 1' Smooth Bore Tubing	50/cs
1056-7-50	Salter Cannula with 7' Smooth Bore Tubing	50/cs
1600-1-50	Salter Cannula with 1' 3-Channel Supply Tubing	50/cs
1600-50	Salter Cannula with 50' 3-Channel Supply Tubing	20/cs
1600Q	Salter Quiet Cannula with 7' 3-Channel Supply Tubing	50/cs
1600TLC	Salter TLC Cannula with Foam Cushions and 7' 3-Channel Supply Tubing	25/cs
1600TLC-1-25	Salter TLC Cannula with Foam Cushions and 1' 3-Channel Supply Tubing	25/cs
1601	Salter Infant Cannula with 1' 3-Channel Tubing	50/cs
1601-7-50	Salter Infant Cannula with 7' 3-Channel Tubing	50/cs
1615	Salter Intermediate Infant Cannula with 7' 3-Channel Supply Tubing	50/cs
1616	Salter Micro Low-Flow Cannula with 7' 3-Channel Supply Tubing	50/cs
1616-4	Salter Micro Low-Flow Cannula with 4' 3-Channel Supply Tubing	50/cs
1616-TLC	Salter Micro TLC Cannula with Foam Cushions and 7' 3-Channel Supply Tubing	25/cs
4824-10	Salter Pediatric Cannula with 7' 3-Channel Supply Tubing	10/cs
4804	Adult Demand Cannula with 4' Supply Tubing	25/cs
4807	Adult Demand Cannula with 7' Supply Tubing	25/cs

## Oxygen Masks

Roscoe Medical's Oxygen Mask is designed for patient comfort with a soft, anatomical shape. Barbed fittings help prevent patient from disconnecting the mask while elastic straps and adjustable aluminum nose clips secure it for a better fit. Standard with crush resistant safety tubing. Latex-free. HCPCS Code: A7015\*.



Item #	Description	Qty
O2MASK-ADULT	Roscoe Adult Medium-Concentration Oxygen Mask with 7' crush-resistant tubing	50/cs
8110	Salter Adult Medium-Concentration Oxygen Mask with 7' three channel supply tubing	50/cs
8030	Salter Adult Elongated High-Concentration Non-rebreather Mask with 7' three channel supply tubing and safety vent	50/cs
8150	Salter Percent O2 Lock Venturi Mask with 7' 3-Channel Supply Tubing	50/cs
0101	Westmed BiFlo Nasal Oxygen Mask with 7' Crush Resistant Supply Tubing	50/cs
1122	Salter Pediatric Oxygen Mask with 7' 3-Channel Supply Tubing	50/cs
8115-0-50	Salter Medium Concentration Oxygen Mask without Supply Tubing	50/cs



## Oxygen Supply Tubing

### Roscoe Medical Oxygen Supply Tubing (Crush-Resistant)

Roscoe Medical's Oxygen Supply Tubing has a three-channel inner lumen design providing uninterrupted flow of oxygen that provides maximum safety for the patient. When uncoiled, tubing lies out straight, minimizing "torque" or pull on the cannula. Our tubing, made without natural rubber latex is available in clear or a distinctive green tint color providing greater visual contrast for improved safety. Universal connectors at each end come standard. Available in a variety of lengths.



TUB-ROS7

Item #	Description	Qty
TUB-ROS7	Three-Channel Crush-Resistant Oxygen Tubing, 7'	50/cs
TUB-ROS25	Three-Channel Crush-Resistant Oxygen Tubing, 25'	25/cs
TUB-ROS50	Three-Channel Crush-Resistant Oxygen Tubing, 50'	20/cs
TUB-ROS25G	Three-Channel Crush-Resistant Oxygen Tubing, 25', Green	25/cs
TUB-ROS50G	Three-Channel Crush-Resistant Oxygen Tubing, 50', Green	20/cs

### Salter Oxygen Tubing

Salter's Three-Channel Safety Supply Tubing has three inside channels that help prevent the tubing from becoming occluded. Once uncoiled, the tubing remains straight to avoid pull on cannula.



2050

Item #	Description	Qty
2014	Three-Channel Safety Tubing, 14'	50/cs
2014G	Three-Channel Safety Tubing, 14', Green	50/cs
2025	Three-Channel Safety Tubing, 25'	25/cs
2025G	Three-Channel Safety Tubing, 25', Green	25/cs
2050	Three-Channel Safety Tubing, 50'	20/cs
2050G	Three-Channel Safety Tubing, 50', Green	20/cs
2001-1-50	Salter 1' 3-Channel Supply Tubing, Clear	50/cs
2002	Salter 7' 3-Channel Supply Tubing, Clear	50/cs
2002G	Salter 7' 3-Channel Supply Tubing, Green	50/cs
2040G	Salter 40' 3-Channel Supply Tubing, Green	20/cs

## Humidifier Adaptor Tubes

Item #	Description	Qty
9996-1	Humidifier Adapter Tubing with Trumpet / Male Thread Grip, 12"	25/cs
SO-676	Humidifier Adapter Tubing with Trumpet / Male Thread Grip, 15"	50/cs
SO1790	Humidifier Adapter Tubing with Trumpet / Male Thread Grip, 21"	50/cs
SHAK-I	Humidifier Adapter with Male Trumpet Thread Grip, 9"	50/cs

## Oxygen Supply Accessories



Item #	Description	Qty
SWV-ROS	Roscoe Medical 360° Rotating Tubing Swivel Connector, Purple	50/cs
WTR-ROS	Roscoe Medical In-Line Water Condensation Trap, Clear	50/cs
0043	Westmed Oxygen Tubing Connector, Barbed, Non-Swivel	50/cs
1016	Salter E-Z Wraps™ Over-The-Ear Foam Tube Cushions, Charcoal, One Pair	50/cs
1220	Salter Tubing Swivel Connector	25/cs
SH75-50	Christmas Tree Adaptor, Tapered And Ribbed, Allows for Secure Tubing Connection	50/pk
1109-2221	Oxysafe Bidirectional Kit	100/cs
1109-2014	2 Inch Connector Tubing for use with Bidirectional Oxysafe	100/cs
1005	Salter Tender Grip Fixation System	25/cs
7001-25	Salter Inline Water Trap with Threaded End Caps	25/cs
SHAK-I	Humidifier Adapter with Male Trumpet Thread Grip, 9"	each
CE1000	Cann-Ease Nasal Moisturizer	each
1225	Salter Swivel Connector with Male/Female Ends	25/cs
1227-2-2-10	Salter Y Connector with Two 2' Bonded Green Tubes	10/cs
1227-7-7-10	Salter Y Connector with Two 7' Bonded Green Tubes	10/cs

## Oxygen Supply Accessories (Cont'd)



SH74



SH749



AG1039642

Item #	Description	Qty
B-F65100	Medline Oxygen Tubing Connector	200/cs
AG1105177	HomeFill Connector Set with Straight Supply Tubing	each
Q51462	Barbed Oxygen Connector	25/pk
SH74	Christmas Tree Adapter, Non-Swivel	50/pk
SH749	Oxygen Y Connector, Barbed	each
AG1039642	Humidifier Elbow Adapter for Everflo	each

## Humidifiers



HUM-ROS6



7900

Item #	Description	Qty
HUM-ROS6	Roscoe Medical Empty Humidifier with 6 Psi Pop-Off Alarm, 350cc,	50/cs
7600	Salter Empty Humidifier with 6 Psi Pop-Off Alarm, 350cc,	50/cs
7100	Salter Empty Humidifier with 3 - 5 Psi Pop-Off Alarm, 350cc,	50/cs
7900-0-25	Salter Empty High Flow Humidifier, 6-15 Lpm with Pop-Off Alarm, 350cc,	25/cs

## Cabinet Filters




CF-1011



CF-1029



CF-1030



CF-1034



CF-1037

Item #	Description	Qty
CF-1000	Cabinet Filter, Roscoe truAIRE-5	20/pk
CF-1011	Cabinet Filter, Invacare MobilAire 5LA, 5LX, 501 and HomeFill	20/pk
CF-1029	Cabinet Filter, DeVilbiss 505 Series	20/pk
CF-1030	Cabinet Filter, Respironics Millennium	20/pk
CF-1033	Cabinet Filter, Airsep NewLife/Elite	20/pk
CF-1034	Cabinet Filter, Invacare Platinum	20/pk
CF-1037	Cabinet Filter, Invacare PerfectO2	20/pk

## Intake and Final Bacteria Filters



Item #	Description	Qty
CIF-O2C5L	Intake Filter, Roscoe truAIRE-5, Rhythm Medical, Direct Supply Concentrators	each
CIF-DEVOEL	Intake Filter, DeVilbiss Extended Life (compare to MC44D-605)	each
CIF-DEVOMF	Intake Filter with Silencer, DeVilbiss Extended Life (compare to 515DZ-605)	each
CIF-EVRFLO	Intake Filter, Respironics EverFlo (compare to 1038831)	each
CIF-IRO	Intake Filter, Invacare Concentrators (compare to 1108078)	each
IBF-RE	Intake Filter, Respironics Millennium (compare to H620)	each
AG1040395	Long Life Filter, Respironics M10 (compare to 1040395)	each
CIF-1127	Pre-Inlet Filter, Respironics Millennium (compare to H618)	20/pk
ILB-316	Final Bacteria Filter, 3/16" Straight Barb	each
ILB-516S	Final Bacteria Filter, 5/16" Straight Barb	each

## Pulse Oximeters

### Roscoe OTC Fingertip Pulse Oximeter

The Roscoe OTC Fingertip Pulse Oximeter allows you to monitor and maintain pulse and oxygen saturation levels on-the-go. Accurate and fast oxygen level and heart rate measuring technology is ideal for sports enthusiasts during strenuous activities or mountain climbing. Recommended for spot-check monitoring only. No prescription required.



Item #	Description	Qty
POX-ROS	Roscoe OTC Fingertip Pulse Oximeter with Lanyard	each

## Analyzers and Sensors

### Salter Analyzers

The Salter Labs Pro2 Check Elite oxygen analyzer will allow you to quickly and easily obtain an accurate measurement of the oxygen purity, Liter Flow and Pressure of stationary and portable oxygen devices. Its ultrasonic technology gives you the accuracy you need without the additional maintenance costs or requiring replacement sensors. It is calibrated out of the box, and is not affected by changing temperatures or altitudes. The Pro2 Check Elite analyzer is equipped with a 5-minute automatic shut-off function to preserve battery life. It includes a carrying case and runs off of a single 9-Volt battery.



Item #	Description	Qty
PRO-2 ELITE-1	Salter PRO2 Check Elite Multifunction Oxygen Analyzer	each

The Maxtec O2+A Oxygen Analyzer gives an accurate oxygen percentage reading in a simple and portable design. It will calibrate to room air or 100% oxygen in one simple touch. Ergonomic design fits comfortably in your hand. Simple one-touch calibration and ON/OFF control for ease of use. Easy-to-read LCD display. Quickly adapts to bore tube applications.



O2A-MAXO2

Item #	Description	Qty
O2A-MAXO2	Maxtec O2+ Oxygen Analyzer	each
O2C-MAXO2	Internal Sensor for Maxtec O2+ A	each

### Equipment Bags and Tags



BAG-15924



BAG-251530R



BAG-5045R



BAG-SETUP



90371

Item #	Description	Qty
BAG-1015	General Purpose Equipment Bag, 10"X15" (1 Mil)	1000/cs
BAG-15924	Nebulizer/Suction Bag, 15"X 9"X24" (1 Mil)	500/cs
BAG-211330R	Concentrator Bag, 21"X13"X30" (1 Mil)	500/rl
BAG-2430R	Equipment Cover, 24"X30" (1.5 Mil)	500/rl
BAG-251530R	Concentrator Bag, 25"X15"X30" (1.5 Mil)	250/rl
BAG-2636R	Equipment Cover, 26"X36" (1.5 Mil)	250/rl
BAG-302035R	Wheelchair/Walker/Commode Bag, 30"X20"X35" (1.5 Mil)	200/rl
BAG-3042R	Equipment Cover, 30"X42" (1.5 Mil)	250/rl
BAG-5045R	Wheelchair Bag, 50"X45" (1 Mil)	250/rl
BAG-SETUP	Patient Set-Up Bag with 4" Gusset Drawstring, 16"X12"	500/cs
90365	Equipment Tag, Dirty (Red)	500/cs
90367	Equipment Tag, Clean (Green)	500/cs
90369	Equipment Tag, Repair (Yellow)	500/cs
90371	Multipack Equipment Tags, 250 each of Clean, Dirty and Repair	750/cs

## Notes

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.



## Notes

[illegible]

## Notes

[illegible]

## This image shows a single page of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins or other markings visible.



6753 Engle Road • Middleburg Hts., OH 44130 Toll Free: (800) 376-7263 • Fax: 440-572-4261

[www.compasshealthbrands.com](http://www.compasshealthbrands.com)

rev\_03\_2023