



THERAPEUTIC SUPPORT SURFACES CATALOG



GROUP I & II Support Surface Qualifiers CMS Coding Eligibility



Medicare Criteria: This information is not intended to be, nor should it be considered, billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare programs, and should consult an attorney or other advisor to discuss specific situations in further detail.

THERAPEUTIC SUPPORT SURFACES

| | |
|--|----|
| Homecare Mattresses | 4 |
| Pressure Redistribution Foam Mattresses | 5 |
| Gel Overlays | 6 |
| Alternating Pressure Pump & Pad Systems | 6 |
| Alternating Pressure and Low Air Loss Mattress Systems | 8 |
| Perimeter Covers | 13 |
| Fall Prevention (Rail Pads & Floor Mats) | 14 |
| Wheelchair Cushions | 15 |



WARNING

Whenever a support surface and bed rails are used as part of any hospital or home care bed system, the potential for bed rail entrapment exists. To reduce the risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment, contact your Roscoe Medical sales representative for a copy of our "Bed Rail Safety Guide" or visit the U.S. FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment".



HPCPS
E0272



Fiber Core Homecare Mattress

- 100% virgin densified polyester fiber core is more economical than foam
- Lighter weight than innerspring mattress, reduces shipping costs
- Breathable polypropylene cover keeps patient cool and comfortable
- Mattress and nylon cover are odor free and non-toxic
- Polypropylene cover is fluid, stain and mold resistant, as well as 100% recyclable
- Meets Federal Fire Code 16 CFR 1633
- One year warranty



| Item # | Description | Size | Wt. Cap. | UOM |
|-----------|--|----------------|----------|------|
| PBHBMATFB | Fiber Core Mattress with Polypropylene Cover | 80" x 35" x 6" | 250 lbs. | each |

HPCPS
E0271



Firm Support Innerspring Mattress



- For use in home care, nursing homes and hospitals where firm support is needed
- Easy-to-clean vinyl cover is durable and fluid-resistant
- High quality 312 bonnell coil provides resilient support base
- 2oz densified polyester fiber topper adds firmness and support
- Heavy-duty resinated pad with a specialized mesh woven back
- Meets 16 CFR 1633 flammability standard
- One year warranty



| Item # | Description | Size | Wt. Cap. | UOM |
|-----------|-----------------------------------|----------------|----------|------|
| PBHBMATIS | Firm Support Innerspring Mattress | 80" x 36" x 6" | 300 lbs. | each |
| PBHBMAT84 | Firm Support Innerspring Mattress | 84" x 36" x 6" | 300 lbs | each |

Bariatric Foam Homecare Mattress



- High density foam core designed specifically for the home care environment
- Durable, long-lasting foam
- Easy-to-clean vinyl cover is durable and fluid-resistant
- Meets Federal Fire Code 16 CFR 1633
- One year warranty



| Item # | Description | Size | Wt. Cap. | UOM |
|-----------|--------------------------------------|----------------|----------|------|
| PBHB42MAT | Bariatric High Density Foam Mattress | 80" x 42" x 7" | 650 lbs. | each |



HCPCS
E0184



SATINFOAM Group 1 Compressed Foam Mattress


GROUP 1

- Premium quality foam homecare mattress
- Easy-to-clean, sheer-resistant nylon cover is durable and fluid-resistant
- Non-slip underside keeps mattress firmly in place
- Super compressed packaging; box footprint is only 14" x 14" x 37"
- Meets 16 CFR 1633 flammability standard
- Limited 12-month warranty



| Item # | Description | Size | Wt. Cap. | UOM |
|-----------------|-----------------------|----------------|----------|------|
| PBSATINFOAM3680 | Group 1 Foam Mattress | 80" x 36" x 6" | 300 lbs. | each |

HCPCS
E0184



Aruba Group 1 Contour Compressed Mattress


GROUP 1

- The sculptured foam surface is designed to provide pressure relief and comfort in three areas: the head and shoulders, the torso and the lower extremities
- The cross-cut design maximizes airflow and provides comfort, support and pressure redistribution
- Deluxe nylon cover is durable, zippered, waterproof and anti-microbial. The nylon top cover reduces friction, heat and contamination. The bottom cover is reinforced vinyl with non-slip protection for safety
- Meets 16 CFR 1633 flammability standard
- Two year warranty on foam, one year on cover



| Item # | Description | Size | Wt. Cap. | UOM |
|-----------|---------------------|----------------|----------|------|
| ARUBA-200 | Aruba Foam Mattress | 80" x 35" x 6" | 350 lbs. | each |

Available
Fall 2023



Aruba-42 Group 1 Bariatric Compressed Foam Mattress

GROUP 1

- Compressed Group 1 Foam to help reduce and prevent pressure sores
- Multi-density foam layers for patient comfort
- Smaller footprint for easier delivery and shipping savings
- Supports those up to 600 lbs.
- Water-resistant zippered cover for easier cleaning
- One year limited warranty

| Item # | Description | Size | Wt. Cap. | UOM |
|----------|-------------------------------|----------------|----------|------|
| ARUBA-42 | Aruba Bariatric Foam Mattress | 80" x 42" x 6" | 600 lbs. | each |



Gel-Lite Gel/Foam Overlay

GROUP 1

- Independent gel tubes prevent migration of gel in high pressure areas
- Raise beds top and bottom without gel migration
- Triple seal polyurethane gel tube construction ensures seal integrity
- High-density foam perimeter on all sides helps prevent patient roll-out
- Conforms to body contours and allows air circulation to reduce heat
- Convuluted foam topper provides additional comfort
- Cover is vapor permeable and moisture resistant; low shear/low friction
- Weight bearing capacity: 250 lbs
- Meets flammability standards TB 117-2013
- One year warranty



| Item # | Description | Size | Wt. Cap. | UOM |
|------------|------------------|----------------|----------|------|
| GO-78343-C | Gel-Lite Overlay | 78" x 34" x 3" | 250 lbs. | each |

Group I | Alternating Pressure Pump and Pad Systems



Meridian SatinAir™ APP System

GROUP 1

- Interlocking bubble cells exchange air every 5 minutes
- Heavy duty 16 gauge PVC
- End flaps secure pad to mattress
- High capacity output 5 liters/minute
- Variable pressure setting allows for maximum comfort
- Adjustable pressure range
- Cycle time: 5 minutes
- Illuminated on/off switch
- Built-in bed mount hooks
- HCPCS Codes: E0181 (pad and pump)
- HCPCS Codes: A4640 (replacement pad)
- One-year limited warranty on pump, 90 days on mattress pad



| Item # | Description | Size | Wt. Cap. | UOM |
|----------|----------------------------------|------------------|----------|-------|
| MER-APP2 | Meridian SatinAir APP System | 78" x 34" x 2.5" | 300 lbs. | each |
| APP-BPD | Basic Replacement Pad | 78" x 34" x 2.5" | 300 lbs. | 10/cs |
| APP-DLXP | Basic Replacement Pad with Flaps | 78" x 34" x 2.5" | 300 lbs. | 10/cs |

A photograph of an elderly man with grey hair, smiling broadly while lying in a hospital bed. He is wearing a light blue hospital gown. A woman with grey hair, wearing a bright pink shirt, stands beside him, also smiling. The background shows medical equipment and a hospital room setting.

CHOOSING THE RIGHT SURFACE

Pressure ulcers are high on the list of "never events" developed by the Center for Medicare and Medicare (CMS) stating that these events should never happen. We carry a broad line of systems with superior quality and exceptional value to meet your patients' individual needs.

HCPCS
E0277



Meridian SatinAir™ 8" APM

GROUP 2

ALTERNATING PRESSURE AND LOW AIR LOSS MATTRESS SYSTEM

- Includes mattress, cover and pump
- Pillow function maintains the air in cells at the head for patient comfort
- Air cells can be easily removed and replaced
- CPR quick deflate at head of mattress
- Quilted nylon zipper cover is low shear, fluid-resistant, and vapor permeable
- 8 liter-per-minute pump in an easy-to-use design
- Two modes of therapy - static or alternating
- Variable pressure setting for maximum comfort
- Support patients up to 350 lbs.
- One-year limited warranty on mattress and pump



| Item # | Description | Size | Wt. Cap. | UOM |
|------------|---------------------------|----------------|----------|------|
| MER-SAPM2 | Meridian SatinAir™ 8" APM | 80" x 36" x 8" | 350 lbs. | each |
| X1SAPM201 | Replacement Mattress | 80" x 36" x 8" | 350 lbs. | each |
| PBSATCOVER | Replacement Cover | 80" x 36" x 8" | N/A | each |

HCPCS
E0277



Meridian SatinAir™ 5" + 3" Foam Base

GROUP 2

5" ALTERNATING PRESSURE AND LOW AIR LOSS MATTRESS SYSTEM WITH 3" FOAM BASE

- Includes mattress, cover and pump
- 3-inch foam base prevents "bottoming out" in the event of a power outage
- Pillow function maintains the air in cells at the head for patient comfort
- Air cells can be easily removed and replaced
- CPR quick deflate at head of mattress
- Quilted nylon zipper cover is low shear, fluid-resistant, and vapor permeable
- 8 liter-per-minute pump in an easy-to-use design
- Two modes of therapy - static or alternating
- Variable pressure setting for maximum comfort
- Pump features both audio and visual alarms
- Support patients up to 350 lbs.
- One-year limited warranty on mattress and pump



| Item # | Description | Size | Wt. Cap. | UOM |
|------------|------------------------|----------------|----------|------|
| MER-SAPM53 | Meridian SatinAir™ 5-3 | 80" x 36" x 8" | 350 lbs. | each |
| X1SAPM5301 | Replacement Mattress | 80" x 36" x 8" | 350 lbs. | each |
| X1SAPM5302 | Replacement Cover | 80" x 36" x 8" | N/A | each |



GROUP 2

Ultra-Care™ 5800

BASIC 8" ALTERNATING PRESSURE/LOW-AIR LOSS MATTRESS SYSTEM

- Includes mattress, cover and pump
- Pillow function maintains the air cells at the head for patient comfort
- Ventilated cell construction for low air loss
- CPR quick deflate at head of mattress
- Air cells made with nylon/PVC
- Quilted cover is water-resistant, vapor-permeable nylon; low shear/low friction
- 8 LPM pump in an easy to use design
- 10 minute cycle times
- Two modes of therapy - static and alternating
- Variable pressure setting allows for maximum comfort
- Support patients up to 350 lbs.
- Two-year limited warranty on mattress and pump



| Item # | Description | Size | Wt. Cap. | UOM |
|-----------|-----------------------------|----------------|----------|------|
| MER-5800 | 8" Mattress with 8 LPM Pump | 80" x 36" x 8" | 350 lbs. | each |
| 5802 | Replacement Mattress | 80" x 36" x 8" | 350 lbs. | each |
| 4802-C-88 | Replacement Cover | 80" x 36" x 8" | N/A | each |

HCPCS
E0277

ULTRACARE



Ultra-Care™ 4800

STANDARD ALTERNATING PRESSURE/LOW AIR LOSS MATTRESS SYSTEM

- Includes mattress, cover and pump
- 8" alternating pressure mattress with low air loss
- Quick disconnect cells for easy cell replacement
- Air cells made with nylon/polyurethane for increased durability and reliability
- CPR quick deflate located at the head of the mattress
- Quilted cover, moisture resistant, vapor permeable nylon; low shear/low friction
- 8 LPM pump in an easy-to-use design
- 10 minute cycle time
- Front panel switch provides two modes of therapy - static or alternating
- Variable pressure setting allows for maximum comfort
- Adjustable pressure range: 20 mmHg – 70 mmHg
- Supports patients up to 350 lbs.
- Two-year limited warranty on mattress and pump



GROUP 2



| Item # | Description | Size | Wt. Cap. | UOM |
|-----------|---|----------------|----------|------|
| 4800 | Ultra-Care 4800 System | 80" x 36" x 8" | 350 lbs. | each |
| 4800E | Ultra-Care 4800 System with 12 LPM Pump | 80" x 36" x 8" | 450 lbs. | each |
| 4802 DC | Replacement Mattress | 80" x 36" x 8" | 450 lbs. | each |
| 4802-C-88 | Replacement Cover | 80" x 36" x 8" | N/A | each |

HCPCS
E0277

ULTRACARE





GROUP 2

HPCPS
E0277



ULTRACARE



Ultra-Care™ 4500

ALTERNATING PRESSURE/LOW-AIR LOSS MATTRESS SYSTEM WITH 3-INCH INTEGRATED FOAM BASE

- System includes mattress, cover and pump
- Features a 5" alternating pressure upper surface combined with an integrated 3" foam base to provide a full eight inches of therapeutic support with bottom-out protection in the event of a power failure
- Separate compartment for foam base inside polyurethane cover prevents soiling
- Quick disconnect cells for easy cell replacement
- Air cells made with nylon/polyurethane for increased durability and reliability
- CPR quick deflate located at the head of the mattress
- Quilted cover, moisture resistant, vapor permeable nylon; low shear/low friction
- 8 LPM pump in an easy-to-use design
- 10 minute cycle time
- Front panel switch provides two modes of therapy - static or alternating
- Variable pressure setting allows for maximum comfort
- Adjustable pressure range: 20 mmHg – 70 mmHg
- Two-year limited warranty on mattress and pump



| Item # | Description | Size | Wt. Cap. | UOM |
|-----------|--|------------------|----------|------|
| MER-4500 | Ultra-Care 4500 Excel | 80" x 35.5" x 8" | 325 lbs. | each |
| 4500E | Ultra-Care Excel System with 12 LPM Pump | 80" x 35.5" x 8" | 425 lbs. | each |
| 4502 | Replacement Mattress | 80" x 36" x 8" | 425 lbs. | each |
| 4502-C-88 | Replacement Cover | 80" x 36" x 8" | N/A | each |

Standard Alternating Pressure Pump



GROUP 2

The Manual Pump for APMs offers both alternating and static functions to provide patients with optimal treatment of pressure ulcers and maximum comfort. It flows at 8 liters-per-minute in a simple, easy-to-use design.



- Front panel switch provides two modes of therapy - static or alternating - for maximum comfort
- LED backlit switches with dial on panel for simple set up and operation
- Includes retractable hooks for bed use
- 8 LPM pump
- Supports patients up to 350 lbs
- Two-year limited warranty



| Item # | Description | UOM |
|--------|-------------------|------|
| 4581 | Standard APM Pump | each |



GROUP 2

**HCPCS
E0277**



ULTRACARE

Ultra-Care™ Xtra

BARIATRIC LOW-AIR LOSS WITH CELL-ON-CELL MATTRESS DESIGN

- Includes mattress, cover and pump
- A combination mattress system providing both alternating pressure and low air loss with a cell-on-cell mattress design for the for patients up to 800 lbs
- Cell-on-cell design prevents the patient from "bottoming out" in power outage
- Quick disconnect cells for easy cell replacement
- 10" air cells made with nylon/polyurethane for increased durability and reliability
- CPR quick deflate located at the head of the mattress
- Quilted cover, moisture resistant, vapor permeable nylon; low shear/low friction
- Pump is a high-flow 12 liter-per-minute pump in an easy-to-use design
- Dual compressor for increased patient comfort and weight capacity
- 10, 15 or 20 minute cycle times
- Two modes of therapy - static or alternating
- Five comfort settings for easy pressure adjustment
- Automatic reset from static to alternating after 45 minutes
- Adjustable pressure range: 20 mmHg – 50 mmHg
- Two-year limited warranty on mattress and pump



| Item # | Description | Size | Wt. Cap. | UOM |
|------------------|--------------------------------------|-----------------|----------|------|
| 4840-42 | Ultra Care Xtra Bariatric APM System | 80" x 42" x 10" | 800 lbs. | each |
| 4840-42-M | Replacement Mattress | 80" x 42" x 10" | 800 lbs. | each |
| PRT-42BARI-COVER | Replacement Cover | 80" x 42" x 10" | N/A | each |

12LPM High Capacity Digital Low Air Loss/Alternating Pressure Pump



GROUP 2



- Upgrade APM system with deluxe high-flow, 12 liters-per-minute air output
- Features dual compressor for increased patient comfort and an added 100 lbs to weight capacity
- Offers both alternating pressure or static operation modes; will automatically reset from static to alternating after 45 minutes
- Includes 10, 15, or 20-minute cycle times and five comfort settings
- Has two ports and is compatible with most Ultra-Care APMs
- Two-year limited warranty



| Item # | Description | UOM |
|--------|------------------------------|------|
| 4253E | Electronic APM Pump – 12 LPM | each |

Magic-Air® Pump

The Quick-Fill Pump is designed to quickly inflate and deflate Low Air Loss Mattresses reducing your patient set-up time by 65%, and allowing you time to attend to your patient's other needs. This quick 20-liter-per-minute flow pump allows your technician to get the mattress up and therapeutic in minutes.



- Light weight, compact and all-steel construction
- Pumps air at a rate of 20-liters-per-minute
- Quieter as compared to standard plastic air pumps
- Weight: 3 lbs.
- Dimensions: 13" length x 5.5" wide height
- 10 foot cord
- Works with Quick Fill Valve 90153 (sold separately)
- One-year limited warranty

| Item # | Description | UOM |
|--------|------------------|------|
| 90501 | Magic-Air® Pump | each |
| 90153 | Quick Fill Valve | each |



Safe-T-Guard™ Cover

PERIMETER COVER TO REDUCE RISK OF FALLS AND ENTRAPMENT

- Reduces the risk of falls and entrapment
- Fits most standard-size hospital mattresses as well as standard eight inch alternating air mattresses
- Secures with three quick-release straps (positioned top, bottom and middle) that connect underneath the mattress
- The quilted top is water resistant and has a covered foam border that rises approximately four inches above the mattress surface and is sloped inward to reduce the risk of roll outs
- Includes cover and 4 foam inserts with separation to allow for patient egress
- Nylon, polyurethane construction to reduce heat build-up
- One-year limited warranty

SGU UNIVERSAL DESIGN



- Secures to mattress with three quick-release straps positioned top, bottom and middle for a snug fit
- Fits most standard hospital bed mattresses 36" W x 80" L x 6" H

SG8 DESIGN (FOR 8" APM SYSTEMS)



- Fits most 8" air mattress systems
- Low shear, nylon zippered cover



| Item # | Description | L x W | Qty. |
|----------|---|-------------|------|
| SGU_FG | Universal Safe-T-Guard Cover with QR Straps | 80" x 35.5" | each |
| SG8_FG | Safe-T-Guard Zipper Cover for 8" Air | 80" x 35.5" | each |
| SGU42_FG | Safe-T-Guard Cover for 42" Mattress | 80" x 42" | each |



LOOKING FOR SOMETHING SPECIFIC?

Your Compass Health Representative is a knowledgeable, trained professional. If you can't find what you need, assistance is just a phone call away. Call **1-800-376-7263** or visit our online catalog at **www.compasshealthbrands.com**.



Bed Rail Pads

- Protects against bumps, bruises
- Vinyl cover wipes clean for easy disinfecting
- Pads do not interfere with the operation of the bed rails
- Attach with hook & loop straps
- One-year limited warranty



| Item # | Description | Size | UOM |
|--------|---------------------------|----------------|------|
| 30038 | Full-Length Bed Rail Pads | 48" x 15" x 1" | pair |
| BR50 | Half-Rail Bed Pads | 31" x 16" x 1" | pair |



High Density Bedside Foam Fall Mat

- Help reduce the possibility of injuries due to falls
- Easy-to-clean vinyl cover
- High density, high compression polyurethane foam
- Flame retardant meets with California Technical Bulletin #117
- Bi-Fold design held together with hook & loop fasteners
- One-year limited warranty



| Item # | Description | Size | UOM |
|--------|----------------------------|----------------|------|
| 30036 | High Density Foam Fall Mat | 66" x 24" x 2" | each |



Vinyl Polymer Bedside Fall Mat with Beveled Edges



- Closed-cell sponge design
- Low profile beveled edge to prevent tripping
- Meets flammability and head impact criteria
- One-year limited warranty
- Proudly made in the USA




| Item # | Description | Size | UOM |
|--------|-------------------------------------|------------------|------|
| FM200 | Bedside Fall Mat with Beveled Edges | 70" x 29" x 5/8" | each |

WHEELCHAIR CUSHIONS

If a patient has a wheelchair for which they meet Medicare coverage criteria, they also qualify for a General Use Seat and Back Cushion (see page 18).



HCPCS
E2601




Meridian General Use Foam Cushion

- Economical seating surface provides comfort and pressure reduction
- Designed to assist in the prevention of pressure ulcers
- Cushion contours the seating area
- Removable, waterproof nylon top with non-slip vinyl bottom
- Color: navy
- Weight bearing: See chart below
- One-year limited warranty
- TAA Compliant options available, ask your Compass Health Representative



| ITEM # | SIZE | Wt. Cap. | UOM | HCPCS |
|-----------|----------------|----------|------|-------|
| WCF-163 | 16" x 16" x 3" | 275 lbs. | each | E2601 |
| WCF-183 | 18" x 16" x 3" | 300 lbs. | each | E2601 |
| WCF-203 | 20" x 16" x 3" | 325 lbs. | each | E2601 |
| WCF-2083 | 20" x 18" x 3" | 350 lbs. | each | E2601 |
| WCBF-223 | 22" x 16" x 3" | 400 lbs. | each | E2602 |
| WCBF-2283 | 22" x 18" x 3" | 425 lbs. | each | E2602 |
| WCBF-2483 | 24" x 18" x 3" | 450 lbs. | each | E2602 |
| GEL-182 | 18" x 16" x 2" | | each | |

HCPCS
E2603



Meridian Optimum Comfort Gel Seat Cushion



- Unique channeled gel bladder prevents gel migration in high-pressure areas
- Gel channels provide relief over the entire cushion area
- Foam base provides additional support
- Inner foam core has a liner to protect foam's integrity
- Weight bearing: See chart below
- Proudly made in the USA
- One-year limited warranty



| ITEM # | SIZE | Wt. Cap. | UOM | HCPCS |
|---------------|-----------|----------|------|-------|
| HDC-9B 16X16 | 16" x 16" | 275 lbs. | each | E2603 |
| HDC-9B 18X16 | 18" x 16" | 300 lbs. | each | E2603 |
| HDC-9B 18X18 | 18" x 18" | 325 lbs. | each | E2603 |
| HDC-9B 20 X16 | 20" x 16" | 350 lbs. | each | E2603 |
| HDC-9B 20X18 | 20" x 18" | 400 lbs. | each | E2603 |
| HDC-9B 22X18 | 22" x 18" | 425 lbs. | each | E2604 |
| HDC-9B 24X18 | 24" x 18" | 450 lbs. | each | E2604 |

HCPCS
E2611




Meridian General Use Back Cushion

- Provides upper back and lumbar support
- Designed to reduce lower back pain and muscle tension
- Removable, waterproof nylon cover
- Stabilization board provides enhanced rigidity and support
- Non-skid backing
- Color: Black
- One-year limited warranty



| ITEM # | SIZE | Wt. Cap. | UOM | HCPCS |
|-----------|-----------|----------|------|-------|
| BKF-1617 | 16" x 17" | NA | each | E2611 |
| BKF-1817 | 18" x 17" | NA | each | E2611 |
| BKF-2019 | 20" x 19" | NA | each | E2611 |
| BKFB-2219 | 22" x 19" | NA | each | E2612 |

HPCPS
E2601



Meridian Optimum Pommel Cushion



- Contoured base and pommel
- Four grades of foam
- Inner core has a liner to protect foam
- Cover is a 4 way stretch polyester/polyurethane
- Moisture repellent & low shear/friction
- Weight bearing: See chart below
- Proudly made in the USA
- One-year limited warranty



| ITEM # | SIZE | Wt. Cap. | UOM | HPCPS |
|-------------|----------------|----------|------|-------|
| HDC-POM1616 | 16" x 16" x 3" | 275 lbs. | each | E2601 |
| HDC-POM1816 | 18" x 16" x 3" | 275 lbs. | each | E2601 |
| HDC-POM2016 | 20" x 16" x 3" | 275 lbs. | each | E2601 |

HPCPS
E2601



Meridian Optimum Coccyx Cushion




- Unique design for relief of pressure on coccyx area
- Weight bearing: See chart below
- One-year limited warranty
- Proudly made in the USA



| ITEM # | SIZE | Wt. Cap. | UOM | HPCPS |
|------------|----------------|----------|------|-------|
| HDC-CX1616 | 16" x 16" x 3" | 225 lbs. | each | E2601 |
| HDC-CX1816 | 18" x 16" x 3" | 225 lbs. | each | E2601 |
| HDC-CX2016 | 20" x 16" x 3" | 225 lbs. | each | E2601 |

HPCPS
E2601



Meridian Optimum Wedge Cushion




- Slopes from 4" at front to 2" at back to reduce risk of fall-outs
- Weight bearing: See chart below
- One-year limited warranty
- Proudly made in the USA



| ITEM # | SIZE | Wt. Cap. | UOM | HPCPS |
|------------|----------------|----------|------|-------|
| HDC-WC1616 | 16" x 16" x 3" | 175 lbs. | each | E2601 |
| HDC-WC1816 | 18" x 16" x 3" | 175 lbs. | each | E2601 |
| HDC-WC2016 | 20" x 16" x 3" | 175 lbs. | each | E2601 |

Medicare Wheelchair Seat Cushion CMS Coding Eligibility



Medicare Criteria: This information is not intended to be, nor should it be considered, billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare programs, and should consult an attorney or other advisor to discuss specific situations in further detail.



6753 Engle Road • Middleburg Hts., OH 44130 Toll Free: (800) 376-7263 • Fax: 440-572-4261

www.compasshealthbrands.com

rev_03_2023