



High Vacuum / Low Flow Aspirator

The **Roscoe Aspirator** is the perfect solution for all your suction needs. It is simple to use and has a durable design to withstand the rigors of homecare use. It produces a maximum vacuum of 675 mmHg and features a high-quality vacuum regulator with anti-vibration setting, an anti-bacteria filter design helps prevent contamination, oil-less, permanent lubrication, and a collection bottle with overflow protection. Designed for easy portability at just 7.7 pounds, its easy to carry and easy to clean.



PRODUCT FEATURES:

- Filter design helps prevent contamination
- High quality vacuum regulator and anti-vibration vacuum gauge setting
- Oil-free, permanent lubrication
- Collection bottle with overflow protection
- Superior performance and high reliability for continuous operation
- Protective gauge guard prevents damage
- Easy carry and cleaning
- HCPCS Code: E0600*

OPTIONAL ACCESSORIES:



PRODUCT SPECIFICATIONS:

Dimensions:	13" L x 6" W x 8.25" H
Weight:	7.7 lbs
Maximum Vacuum:	-675 +/- 3% mm Hg
Max Air Flow:	-15-17 LPM
Power Cord Length:	6'
Electrical Requirements:	110-120V/220-240V (50-60Hz)
Pressure Display:	Analog Gauge
Filter Pore Size:	0.45 microns
Collection Capacity:	800 cc
Operating Temperature:	32° - 104°F
Operating Humidity:	95% RH
Warranty:	Two Year Limited



Item #	Description	Qty
ASP-ROS	Roscoe High Vacuum / Low Flow Aspirator	each
90386	Suction Kit (Includes elbow, filter, collection jar, 6' tubing and two 7" tubes)	each
Replacement Parts:		
50038	800 cc Disposable Collection Jar	12/cs
ILB-ROS	Filter, 3/8" to 1/2" barbed	each
ROS-610E	Elbow	20/cs

* This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional. Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.