



Bariatric Bath Bench

The **ProBasics Bariatric Bath Bench** is designed for individuals who require seating support and stability when bathing or showering. The bench features a special cross brace design with a heavy duty blow-molded textured seat that can support patient weights up to 550 pounds. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Comfortable, blow-molded seat with integrated handles for additional support
- Dual cross brace design reinforced with aircraft-style rivets for extra stability
- High-strength aluminum frame to prevent corrosion
- Textured contoured surface with drainage holes to reduce slipping and increase comfort
- Wide-based tips help increase stability
- Legs are height adjustable in one-inch increments with heavy duty height adjustment pins to accommodate larger individuals
- Quick tool-free assembly
- HCPCS Code: E0245*

SPECIFICATIONS

Overall Footprint:	21"W x 17.5"D
Seat Dimensions:	20"W x 12"D
Adjustable Seat Height:	16.5" – 20.5"
Weight Capacity:	550 lb
Carton Dimensions:	18.5" x 22.5" x 17.5"
Shipping Weight:	12 lb
Warranty:	Limited Lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSBC	Bariatric Bath Bench	815067071944	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.