



Infrared No-Touch Forehead Thermometer

The Sejoy Infrared No-Touch Thermometer delivers instant, single second temperature reads on both foreheads and objects alike. With a fever alarm and memory function to store up to 10 temperature readings, the Infrared Thermometer will be your go-to for quick and reliable temperature measurements. Dual scale readings in Fahrenheit and Celsius, low battery indicator and auto shut-off features also included. Retail package comes complete with Infrared Thermometer, 2 AAA Alkaline Batteries and User Manual.

FEATURES

- Fast 1 second readings; beeps when complete
- Memory of up to 10 temperature readings
- Optional backlight for easy reading
- Results indicator lets you know when reading is ready
- Non-touch measurement for forehead or objects such as baby bottles
- Dual scale modes – Fahrenheit or Celsius
- Low battery indicator
- Fever alarm
- Auto shut-off
- Includes 3-Color Backlight display for normal, moderate and high temperature readings

SPECIFICATIONS

Product Dimensions:	6" L x 1.6" W x 1.7" D
Retail Box Dimensions:	6.93" L x 3.66" W x 2.05" D
Product Weight:	2.96 oz (with batteries)
Measure Time:	Approximately 1 second
Measuring Range:	Forehead Mode: 93.2°F ~ 109.4°F (34°C ~ 43°C)
	Object Mode: 32°F ~ 212°F (0°C ~ 100°C)
Battery Size:	AAA Battery (2)
Battery Life:	Approximately 1 year/6000 readings
Expected Service Life:	3 Years
Warranty:	One year limited

Item #	Description	GTIN	UOM
DET-306L	Infrared No-Touch Forehead Thermometer, with 3-Color Backlight	06970392215631	each