



## K4 Lite Wheelchair

### HIGH STRENGTH LIGHTWEIGHT WHEELCHAIR

The **ProBasics K4 Lite Wheelchair** is a high-strength lightweight extremely versatile wheelchair. With an included seat and upholstery extension combined with all of its other unique features that come standard, the ProBasics K4 Lite easily adapts to any patient's individual needs. It offers superior value and functionality for those looking for a high-performance mobility product at an affordable price.



Flip-back armrests are height adjustable



2" upholstery seat extension



Adjustable anti-tippers with push-pin flip-up wheels



Quick release rear axle

### FEATURES

- 3-Position quick-release axles allow for tool-free removal of rear wheels for fast hemi-height adjustment
- Comes included with a seat rail and upholstery extension that easily adjusts seat depth from 16" to 18" without any additional upholstery
- Flip-back padded arm rests are height adjustable from 10" to 14"
- Adjustable height wrap-around back upholstery for maximum comfort
- Comfortable and stylish nylon upholstery
- Polyurethane tires and composite casters for superior performance
- Includes tool-free adjustable front riggings
- Adjustable anti-tippers included
- Easy push-to-lock wheel locks
- HCPCS\* Code: K0004

K0195 or E0990 (elevating legrests only)

E0971 (anti-tippers)

K0040 (angle adjustment footplates)

E2201 (non-standard seat frame - 20" chair only)

### SPECIFICATIONS

Seat Width:	16"   18"   20"
Seat Depth:	16" – 18"
Seat-to-Floor Height <sup>1</sup> :	Adult 19.5"   Hemi 17.5"
Back Height:	16" – 17.5"
Back Style:	Height Adjustable
Overall Width w/Arms:	24.5"   26.5"   28.5"
Overall Depth:	32" (without riggings)
Weight Capacity:	300 lb
Product Weight <sup>2</sup> :	32 – 36 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	49 – 52 lbs (Includes Legrests)
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC4L1616DS	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067078547	each
WC4L1816DS	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067078554	each
WC4L2016DS	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067078561	each
WC4L1616DE	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067078516	each
WC4L1816DE	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067078523	each
WC4L2016DE	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067078530	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

<sup>1</sup>NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

<sup>2</sup>NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.