




Aluminum Transport Chair

The **ProBasics Aluminum Transport Chair** is a compact, lightweight design with an easy-to-maintain finish. The flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport. Comes complete with seat belt, tool-free swing-away composite footrests and heel loops. Supports patient weights up to 300 pounds.




Seat belt



Handy back-release hinge



Flared handles back to fold down flat



Push-to-lock rear wheel locks

FEATURES

- Lightweight aluminum construction (only 20 lbs)
- Flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport
- Stylish aluminum frame with permanent rear axle for strength and durability
- Padded full-length permanent armrests
- Padded nylon upholstery is comfortable, lightweight and easy to clean
- Seat belt for added safety
- Nylon seat straps allow for ease of folding and lifting
- Push-to-lock rear wheel locks
- Composite 8" wheels with solid rubber tires are lightweight and maintenance free
- Tool-free, swing-away composite footrests with heel loops included
- HCPCS Code: E1038*

COLORS

Blue
Finish



Burgundy
Finish



Pink
Finish



SPECIFICATIONS

Overall Width:	22" open 9" closed
Width b/w Grips:	19.75
Overall Height:	37"
Overall Depth:	38.75" with riggings 26.75" without riggings
Seat-to-Floor Height:	19"
Back Height:	16"
Seat Dimensions:	19" (W) x 16" (D)
Weight Capacity:	300 lb
Product Weight:	20 lb
Carton Dimensions:	23" x 11" x 31"
Shipping Weight:	26 lb
Warranty:	Limited Lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
TCA1916BL	ProBasics Aluminum Transport Chair	Blue	815067070503	1/cs
TCA1916BG	ProBasics Aluminum Transport Chair	Burgundy	815067070497	1/cs
TCA1916PK	ProBasics Aluminum Transport Chair	Pink	815067070510	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.