




## Three-Position Recliner


The **PROBASICS®** three-position recliner offers both comfort and durability for day room or patient room use. Featuring full recline, elevated footrest with elevated and full upright positions, this recliner provides maximum flexibility for users. Need to move the recliner? The convenient push bar on the back and four rolling casters let you glide the chair where it's needed! For user safety, the two rear steering casters can be locked to secure the chair in place.



### FEATURES

- Hide-away footrest folds completely below chair to make rising to a standing position easier
- Vinyl seat and back padding to help reduce bottoming out
- Adjustable height settings
- Rust resistant rear locking caster forks
- Easy-clean, blow-molded side panels can be fast with no tools, just press the pin
- 5 position, locking molded trays easily slide ON and OFF with built-in side storage when not in use
- Convenient overlap tray with high lip completes the chair and the tabs let you move the tray freely
- Articulating headrest helps the patient stay within the social "field of vision"
- Flame retardant upholstery – meets California Technical Bulletin 117
- Prop 65 and CARB marked for sale in California

### AVAILABLE COLORS



### SPECIFICATIONS

Overall Height:	47"
Overall Width:	24"
Overall Depth Upright:	40"
Overall Depth Recline:	76"
Seat Width:	19"
Seat Depth:	19.5"
Seat-to-Back Height:	32"
Seat-to-Floor Height:	21.5"
Seat Angle Range	13°-38°
Back Angle Range	28°-54°
Tray Size	25.5" L x 15" D
Weight Capacity:	250 lbs
Product Weight	73 lbs
Shipping Weight:	86 lbs
HCPCS*	E1031

Item #	Description	UPC	UOM
PB6074BL	ProBasics Three-Position Recliner, Blue	815067072903	1/cs
PB6074RW	ProBasics Three-Position Recliner, Rose	815067072910	1/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.