



# Deluxe Height Adjustable Aluminum Rollator


The stylish **ProBasics Deluxe Height Adjustable Aluminum Rollator** features a unique design allowing the user to adjust the seat and handle height to their specifications for maximum comfort and universal fit (adjusts from junior to adult size with one rollator). This lightweight rollator feature easy-to-operate ergonomic hand brakes and 8-inch wheels for great maneuverability and stability. A seamless padded flip-up seat with zipper compartment and large nylon storage pouch keep belongings private and protected.



Water-Resistant Nylon Pouch



Easy-to-Operate Ergonomic Hand Brakes



8-inch Wheels

## FEATURES

- Ergonomic handles with loop-style brakes are adjustable in one-inch increments from 30" - 38" with self-threading knob
- Height adjustment on frame allows seat height to be adjusted in one-inch increments from 18" to 22" (Figure B). Comes with a padded flip-up seat with zipper pouch for storage of personal items and added security.
- Additional roomy and convenient storage pouch attaches to frame easily and keeps larger belongings protected.
- Removable padded curved backrest is equipped with a hinge so it can flip up for easy storage.
- 8-inch non-marring black wheels with soft grip tires glide on almost any surface.
- Lightweight aluminum frame is durable for years of worry-free use.
- Folds up easily for easy transport and storage.
- First-Time Tool-Free Assembly Required
- New and improved carton size for shipping and costing.
- HCPCS\* Code: E0143, E0156 (seat)

## COLORS



## SPECIFICATIONS

Overall Width:	23" (W)
Seat Dimensions:	14" (W) x 12" (D)
Adjustable Seat Height:	18" - 22"
Patient Height Range:	4'11" - 6'4"
Weight Capacity:	300 lb
Adjustable Handle Height:	30" - 38"
Wheel Size:	8"
Product Weight:	14 lb
Carton Dimensions:	21.5" x 7.8" x 20"
Shipping Weight:	19 lb
Warranty:	Limited lifetime on frame

Item #	Description	Color	UPC	UOM
RLAA8BG	ProBasics Adjustable Rollator with 8-inch Wheels	Burgundy	815067072897	1/cs
RLAA8BL	ProBasics Adjustable Rollator with 8-inch Wheels	Blue	815067072880	1/cs
RLAA8BK	ProBasics Adjustable Rollator with 8-inch Wheels	Black	815067072866	1/cs
RLAA8PK	ProBasics Adjustable Rollator with 8-inch Wheels	Pink	815067072873	1/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.