


HCPCS  
E2601



## Meridian General Use Foam Cushion

- Economical seating surface provides comfort and pressure reduction
- Designed to assist in the prevention of pressure ulcers
- Cushion contours the seating area
- Removable, waterproof nylon top with non-slip vinyl bottom
- Color: navy
- Weight bearing: See chart below
- One-year limited warranty
- TAA Compliant options available, ask your  
Compass Health Representative



ITEM #	SIZE	Wt. Cap.	UOM	HCPCS
WCF-163	16" x 16" x 3"	275 lbs.	each	E2601
WCF-183	18" x 16" x 3"	300 lbs.	each	E2601
WCF-203	20" x 16" x 3"	325 lbs.	each	E2601
WCF-2083	20" x 18" x 3"	350 lbs.	each	E2601
WCBF-223	22" x 16" x 3"	400 lbs.	each	E2602
WCBF-2283	22" x 18" x 3"	425 lbs.	each	E2602
WCBF-2483	24" x 18" x 3"	450 lbs.	each	E2602